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Reference	miRNA	Animal	Delivery	Doses	Volume	Injection Times	Detection
Min-Jie Ni, et al. PLOS ONE 2011	mi1-HongES2	7-10 days SD rat	cauda epididymis injection	4nmol per side	N/A	N/A	3 days after injection
Dong Li, et al. JBC 2011	miR-99a	nude mice	tumor mass injection	10nmol	0.1 ml per injection	every three days for two weeks	2-4 weeks after injection
Jin Hou, et al. Cancer Cell 2011	miR-199a/b-3p	nude mice	tumor mass injection	10nmol	0.1 ml per injection	once every three days for two weeks	2-4 weeks after injection
Kutzfeldt J, et al. Nature 2005	miR-122 miR-192	6-week-old mice	tail-vein injections	80 mg per kg body weight	0.2 mL per injection	1 to 3 consecutive days	24 h after the last injection
Jean Krutzfeldt, et al. Nucleic Acids Research 2007	miR-16	6-week-old mice	Intracerebral injection	~800 ng per mouse	0.2 ml per injection	a single injection	72 h after injection
Alessandra Care, et al. Nature Medicine 2007	miR-133	8-week-old mice	minipump implantation	80 mg per kg body weight	N/A	a single infusion	1 month after minipump implantation
Laura Fontana, et al. PLOS ONE 2008	miR-17-5p	6-week-old nude mice	intratumorally injection	0.2 ng per mouse	0.1 ml	three times per week for two weeks	2 weeks after the first injection
Thomas Thurn, et al. Nature 2008	miR- 21	10-12 weeks old mice	jugular vein catheter	80 mg per kg body weight	N/A	injected daily for three days	3 weeks after the first injection
Neri Mercatelli, et al. PLOS ONE 2008	miR-221 miR-222	6-week-old mice	intratumorally injection	1 ug of each antagomir	40 µl	at day 0, 5 and 9, for a injections per tumor	N/A
Ashutosh Dharap, et al. J Cereb Blood Flow Metab 2009	miR-145	rat	minipump implantation	totally 960 pmol each rat	48 µl	N/A	N/A
L.H. Xie H, et al. J Clin Invest 2009	miR-2861	6-week-old female mice	tail vein injection	80 mg/kg body weight	0.2 ml per injection	days 1 to 3 for 3 consecutive weeks	4 days, 3 weeks, and 6 weeks after the last injection
Bonaure A, et al. Science 2009	miR-92	mice	intravenous injection	8 mg per kg of body weight	0.2 ml per injection	at days 0, 2, 4, 7, 9 of the left coronary artery	14 days
Jie Ding, et al. Nature cell biology 2010	miR-151-5p	nude mice	transplantation	2x10 ⁶ cells transfected by antagomir	N/A	N/A	5 weeks after injection
Linfu Liang, et al. Hepatology 2010	miR-125b	6-8 weeks old nude mice	subcutaneously injection	2x10 ⁶ cells transfected by antagomir	N/A	N/A	4 weeks after injection
Mingfeng Shu, et al. Molecular Cancer 2011	miR-335	4-week-old nude mice	intratumorally injection	4 mg/ml	50 µl	every two days for 2 weeks	2 weeks after the first injection