



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

INSTRUCTION MANUAL

ZR Genomic DNA I Kit™

Catalog Nos. **D3004 & D3005**

Highlights

- Easy purification of high quality DNA from whole blood, plasma, serum, body fluids, buffy coat, lymphocytes, tissue, swabs or cultured cells in less than 20 minutes using innovative *ZymoBead™* silica-bead technology.
- Compatible with commonly used anticoagulants (i.e., EDTA, heparin, citrate).
- Excludes the use of Proteinase K and organic denaturants.
- Isolated DNA is ideal for PCR, endonuclease digestion, bisulfite conversion/methylation detection, sequencing, genotyping, etc.

Contents

Product Contents & Specifications.....	1
Product Description	2
Protocols	
Whole Blood, Serum, and Plasma Samples..	3
Buccal Cell Samples.....	4
Solid Tissue Samples.....	5
Cell Monolayer Samples.....	5
Cell Suspensions and Proteinase K Digested Samples.....	7
Troubleshooting.....	8
Ordering Information.....	8
List of DNA-Related Products	9

Satisfaction of all Zymo Research products is guaranteed. If you are not satisfied with this product please call 1-888-882-9682.

Product Contents

ZR Genomic DNA I Kit™	D3004 (100 purifications)	D3005 (400 purifications)	Storage Temperature
Genomic Lysis Buffer	50 ml	2 x 100 ml	Room Temp.
DNA Pre-Wash Buffer	30 ml	2 x 50 ml	Room Temp.
g-DNA Wash Buffer	50 ml	2 x 100 ml	Room Temp.
DNA Elution Buffer	4 ml	16 ml	Room Temp.
ZymoBeads™	1 ml	4 x 1ml	Room Temp.
Instruction Manual	1	1	-

Note - Integrity of kit components is guaranteed for up to one year from date of purchase. Reagents are routinely tested on a lot-to-lot basis to ensure they provide maximal performance and reliability.

Specifications

- **Sample Sources** – Whole blood, plasma, or serum from humans, mice, rats, etc. Also, tissue, cells from culture, buccal cells, as well as a variety of biological liquids are effectively processed using this kit.
- **DNA Purity** – High-quality DNA is eluted with **DNA Elution Buffer** or water. DNA is well suited for PCR and other downstream applications. Typical absorption indices are $A_{260}/A_{280} > 1.7$
- **DNA Size Limits** – Capable of recovering DNA up to and above 40 kb. In most instances, mitochondrial DNA and viral DNA (if present) will also be recovered.
- **ZymoBead™ Binding Capacity** – ~5 µg DNA per 10 µl ZymoBead™ suspension.
- **DNA Recovery** – Typically, DNA is eluted into 35 µl **DNA Elution Buffer** or water for the standard procedure. Human whole blood will yield 3 - 7 µg DNA per 100 µl blood sampled. Mammalian tissues yield: 1 - 3 µg DNA per mg skeletal, heart, and brain tissues and 3 - 5 µg DNA per mg liver, kidney and lung tissues.
- **Equipment** – microcentrifuge, vortex


Note - ™ Trademarks of Zymo Research Corporation. This product is for research use only and should only be used by trained professionals. Some reagents included with this kit are irritants. Wear protective gloves and eye protection. Follow the safety guidelines and rules enacted by your research institution or facility.

ZYMO RESEARCH CORP.

Toll Free: 1-888-882-9682 • Fax: 1-714-288-9643 • Web: www.zymoresearch.com • E-mail: info@zymoresearch.com


Product Description

The **ZR Genomic DNA I Kit™** is a simple procedure for the rapid isolation of total DNA (e.g., genomic, mitochondrial, viral) from a variety of biological sample sources. This product has been optimized for maximal recovery of ultra-pure DNA without RNA contamination and is compatible with whole blood (fresh or stored), serum, plasma, buffy coat, solid tissue, bone marrow and buccal cells, cells from culture, and many biological liquid samples. For processing, simply add the specially formulated **Genomic Lysis Buffer** to a sample in a 1.5 ml tube, add *ZymoBeads™*, vortex, then centrifuge. There is no need for organic denaturants or Proteinase K digestion because of the unique chemistries featured in the kit that yield high-quality, purified DNA in just minutes (see below). PCR inhibitors are effectively removed during the purification process. DNA purified using the **ZR Genomic DNA I Kit™** is suitable for PCR, nucleotide blotting, DNA sequencing, endonuclease digestion, bisulfite conversion/methylation analysis, and other downstream applications.




Ultra-pure DNA is ideal for...

- PCR
- Endonuclease Digestion
- Genotyping
- Bisulfite Conversion & Methylation Analysis



DNA isolation from different sources using the **ZR Genomic DNA I Kit™**. Purifications were performed in duplicate for each sample and an equal volume of eluted DNA was loaded into each lane of a 0.8% (w/v) TAE/agarose/ethidium bromide gel. M is a 1 kb DNA ladder.



Endonuclease digestion of DNA purified with the **ZR Genomic DNA I Kit™**. M is a 1 kb DNA ladder.

The **ZR Genomic DNA II Kit™** (D3006, D3007, D3024, D3025) provides a spin column alternative for isolation of up to 25 µg DNA/column.

Zymo Research offers the **EZ DNA Methylation™** (D5001, D5002, D5003, D5004), **EZ DNA Methylation-Gold™** (D5005, D5006, D5007, D5008) and **EZ DNA Methylation-Direct™** (D5020, D5021, D5022, D5023) Kits for rapid, precise DNA methylation detection.

For Technical Assistance, please contact those at Zymo Research's Technical Department at 1-888-882-9682 or E-mail to tech@zymoresearch.com.

The pellet should be completely resuspended to ensure no impurities are transferred to the next step.

Elution of DNA from the ZymoBeads™ is dependent on pH and temperature. If water is used, ensure the pH is > 6.0. Also, the total yield may be improved by eluting the DNA with Elution Buffer or water pre-equilibrated to 60-70°C.

PROTOCOLS

WHOLE BLOOD, SERUM, AND PLASMA SAMPLES

The following is for the purification of DNA from 50 µl whole blood, serum or plasma; however, the volume can be adjusted depending on your requirements. Fresh, frozen, or preserved blood (in EDTA, citrate, or heparin) can be used. If whole blood cannot be processed immediately, the sample can be “stabilized” for later processing (as noted below) although the immediate processing of blood samples is recommended.

1. Ensure the **ZymoBeads™** are fully resuspended by vortexing. In a 1.5 ml tube, add 200 µl of **Genomic Lysis Buffer** to 50 µl of blood, then add 10 µl **ZymoBeads™**. Mix by inversion, and then incubate at room temperature for 5 minutes. Centrifuge the tube at 1,500 x *g* for 1 minute. Carefully remove the supernatant without disturbing the bead pellet.

Note: For multiple sample processes, make a **Genomic Lysis Buffer-ZymoBeads™** mixture. For example, mix 40 µl **ZymoBeads™** to 800 µl of **Genomic Lysis Buffer**. Add 210 µl of this mixture to each 50 µl blood sample.

2. Add 200 µl of **Genomic Lysis Buffer** to the **ZymoBeads™**. Resuspend pellet by pipetting up and down. Centrifuge at 1,500 x *g* for 1 minute. Discard the supernatant.
3. Add 200 µl of **DNA Pre-Wash Buffer** to the **ZymoBeads™**. Resuspend the pellet and then centrifuge at 1,500 x *g* for 1 minute. Discard the supernatant.
4. Add 500 µl of **g-DNA Wash Buffer** to the **ZymoBeads™**. Resuspend the pellet and then centrifuge at 1,500 x *g* for 1 minute. Discard the supernatant. Recentrifuge briefly and remove any residual wash buffer.
5. Add ≥ 35 µl of **Elution Buffer** or water, resuspend pellet by pipetting up and down, and then centrifuge at 10,000 x *g* for 1 minute.
6. Collect the supernatant. The supernatant contains purified DNA that can be used immediately or stored (e.g., -20°C) for later use.

Delayed Processing of Blood Samples: The immediate processing of blood with this kit is recommended. However, if blood cannot be processed immediately, samples can be “stabilized” in **Genomic Lysis Buffer** for processing at a later time. This can be achieved by adding four volumes of **Genomic Lysis Buffer** to each volume of whole blood (4:1). Blood samples mixed with **Genomic Lysis Buffer** can be stored at room temperature for 7 days, 0 - 4°C for up to 15 days, -20°C for 2 months, or < -70°C for many years. Samples stored at ≤ 4°C should reach room temperature prior to the addition of **ZymoBeads™** as given in the standard protocol (above).

BUCCAL CELL SAMPLES

Buccal cells can be isolated using a rinse- or swab-based isolation method.

- A. **Rinse Method:** Vigorously rinse 10 - 20 ml of saline solution or mouthwash orally for 30 seconds. The more vigorous the rinsing action, the more cells that will be recovered. Spit the solution into a 50 ml tube and pellet the cells at 1,500 x *g* for 5 minutes. Discard the supernatant without disturbing the cell pellet. Add 500 µl of **Genomic Lysis Buffer** to the pellet then vortex for 4 - 6 seconds.
- B. **Swab Isolation Method:** Thoroughly rinse mouth out before isolating cells. Brush the inside of the cheek with a buccal swab for 15 seconds (approximately 20 brushes), making sure to cover the entire area of the inner cheek. Rinse the brush into a microcentrifuge tube using 500 µl of **Genomic Lysis Buffer** then vortex for 4 - 6 seconds.
1. Ensure the **ZymoBeads™** are fully resuspended by vortexing. Add 10 µl **ZymoBeads™** to the **Genomic Lysis Buffer**/cell mixture in A or B (above). Mix by inversion, and then incubate at room temperature for 5 minutes. Centrifuge the tube at 1,500 x *g* for 1 minute. Discard the supernatant.
 2. Add 200 µl of **Genomic Lysis Buffer** to the **ZymoBeads™**. Resuspend pellet by pipetting up and down. Centrifuge at 1,500 x *g* for 1 minute. Discard the supernatant.
 3. Add 200 µl of **DNA Pre-Wash Buffer** to the **ZymoBeads™**. Resuspend the pellet and then centrifuge at 1,500 x *g* for 1 minute. Discard the supernatant.
 4. Add 500 µl of **g-DNA Wash Buffer** to the **ZymoBeads™**. Resuspend the pellet and then centrifuge at 1,500 x *g* for 1 minute. Discard the supernatant. Recentrifuge briefly and remove any residual wash buffer.
 5. Add ≥ 35 µl of **Elution Buffer** or water, resuspend pellet by pipetting up and down, and then centrifuge at 10,000 x *g* for 1 minute.
 6. Collect the supernatant. The supernatant contains purified DNA that can be used immediately or stored (e.g., -20°C) for later use.

Elution of DNA from the ZymoBeads™ is dependent on pH and temperature. If water is used, ensure the pH is > 6.0. Also, the total yield may be improved by eluting the DNA with Elution Buffer or water pre-equilibrated to 60-70°C.

Soft tissue samples are readily homogenized using our **Squisher™-Single**, **Squisher™-8**, and **Squisher™-96** products.

Typical yields are: 1 - 3 µg DNA per mg skeletal, heart, and brain tissues and 3 - 5 µg per mg liver, kidney, and lung tissues.

Elution of DNA from the ZymoBeads™ is dependent on pH and temperature. If water is used, ensure the pH is > 6.0. Also, the total yield may be improved by eluting the DNA with Elution Buffer or water pre-equilibrated to 60-70°C.

SOLID TISSUE SAMPLES

Note: For Proteinase K digested materials (e.g., tailsnips) follow the protocol for **CELL SUSPENSIONS AND PROTEINASE K DIGESTED SAMPLES** (pg. 7). Otherwise...mechanically homogenize up to 5 mg of fresh or frozen tissue in 500 µl of **Genomic Lysis Buffer**. Increase the reagents proportionally if more than 5 mg of solid tissue is used.

1. Centrifuge the lysate at top speed (10,000 x g) for 5 minutes. Making sure not to disturb the pelleted debris, transfer the supernatant to a clean microcentrifuge tube.
2. Ensure the **ZymoBeads™** are fully resuspended by vortexing. Add 10 µl **ZymoBeads™** to the supernatant. Mix by inversion, and then incubate at room temperature for 5 minutes. Centrifuge the tube at 1,500 x g for 1 minute. Discard the supernatant.
3. Add 200 µl of **Genomic Lysis Buffer** to the **ZymoBeads™**. Resuspend pellet by pipetting up and down. Centrifuge at 1,500 x g for 1 minute. Discard the supernatant.
4. Add 200 µl of **DNA Pre-Wash Buffer** to the **ZymoBeads™**. Resuspend the pellet and then centrifuge at 1,500 x g for 1 minute. Discard the supernatant.
5. Add 500 µl of **g-DNA Wash Buffer** to the **ZymoBeads™**. Resuspend the pellet and then centrifuge at 1,500 x g for 1 minute. Discard the supernatant. Recentrifuge briefly and remove any residual wash buffer.
6. Add ≥ 35 µl of **Elution Buffer** or water, resuspend pellet by pipetting up and down, and then centrifuge at 10,000 x g for 1 minute.
7. Collect the supernatant. The supernatant contains purified DNA that can be used immediately or stored (e.g., -20°C) for later use.

CELL MONOLAYER SAMPLES

The following procedure is designed for up to 1.0×10^6 monolayer cells (roughly equal to one well of a 6-well plate or ½ of a T25 flask). Although cell types and culture conditions may vary, the protocol will work with high-density growth cells (e.g., HeLa cells) as well as with low-density growth cells (e.g., neuronal cells). The procedure may be scaled up or down for increases or decreases in the amounts of monolayer cells sampled (see the **Guidelines for Monolayer Cell DNA Isolation** below).

1. Trypsinize or manually scrape adherent cells from the growth surface of a culture flask or plate. Centrifuge the cell suspension at approximately 500 x g for 5 minutes. Remove the supernatant and add 500 µl of **Genomic Lysis Buffer** directly to the pellet. Resuspend pellet by vortexing 4 - 6 seconds.

2. Ensure the **ZymoBeads™** are fully resuspended by vortexing. Add 10 µl **ZymoBeads™** to the supernatant. Mix by inversion, and then incubate at room temperature for 5 minutes. Centrifuge the tube at 1,500 x g for 1 minute. Discard the supernatant.
3. Add 200 µl of **Genomic Lysis Buffer** to the **ZymoBeads™**. Resuspend pellet by pipetting up and down. Centrifuge at 1,500 x g for 1 minute. Discard the supernatant.
4. Add 200 µl of **DNA Pre-Wash Buffer** to the **ZymoBeads™**. Resuspend the pellet and then centrifuge at 1,500 x g for 1 minute. Discard the supernatant.
5. Add 500 µl of **g-DNA Wash Buffer** to the **ZymoBeads™**. Resuspend the pellet and then centrifuge at 1,500 x g for 1 minute. Discard the supernatant. Recentrifuge briefly and remove any residual wash buffer.
6. Add ≥ 35 µl of **Elution Buffer** or water, resuspend pellet by pipetting up and down, and then centrifuge at 10,000 x g for 1 minute.
7. Collect the supernatant. The supernatant contains purified DNA that can be used immediately or stored (e.g., -20°C) for later use.

Elution of DNA from the ZymoBeads™ is dependent on pH and temperature. If water is used, ensure the pH is > 6.0. Also, the total yield may be improved by eluting the DNA with Elution Buffer or water pre-equilibrated to 60-70°C.

Guidelines for Monolayer Cell DNA Isolation: The above procedure is designed for the processing of 0.1 - 1.0 x 10⁶ cells. However, cell numbers (growth densities) can vary between different cell types. Table 1 (below) provides an approximation of what can be recovered from different culture containers for high-density growth cells like CV1 and HeLa cells. The table provides a guideline when adjusting reagent volumes used in the protocol. If the cell number is less than 1.0 x 10⁵, use half the volume of **Genomic Lysis Buffer** given in the standard protocol (above). If the cell number is greater than 1.0 x 10⁶, use double the volume of **Genomic Lysis Buffer** and **ZymoBeads™** indicated in the protocol.

Table 1: Culture Plate/Flask Growth Area (cm²) and Cell Number

Culture Container	Well /Flask Surface Area	Cell Number
96-well plate	0.32 - 0.6 cm ²	4 - 5 x 10 ⁴
24-well plate	2 cm ²	1 - 3 x 10 ⁵
12-well plate	4 cm ²	4 - 5 x 10 ⁵
6 -well plate	9.5 cm ²	0.5 - 1 x 10 ⁶
T25 Culture Flask	25 cm ²	2 - 3 x 10 ⁶
T75 Culture Flask	75 cm ²	0.6 - 1 x 10 ⁷
T175 Culture Flask	175 cm ²	2 - 3 x 10 ⁷

CELL SUSPENSIONS AND PROTEINASE K DIGESTED SAMPLES

Cells should be processed directly from biological fluids or from suspension in PBS, TE, or compatible buffers.

If the cell number is $>1.0 \times 10^6$, increase the volume of ZymoBeads™ added to the sample. (e.g., 20 μ l ZymoBeads™ for 2.0×10^6 cells in 200 μ l PBS)

Elution of DNA from the ZymoBeads™ is dependent on pH and temperature. If water is used, ensure the pH is > 6.0 . Also, the total yield may be improved by eluting the DNA with Elution Buffer or water pre-equilibrated to 60-70°C.

Typical yields from Proteinase K digested tissues are: 1-3 μ g DNA per mg skeletal, heart, and brain tissues and 3-5 μ g per mg liver, kidney, and lung tissues.

The following protocol is designed for up to 200 μ l of biological liquid sample including CSF, buffy coat, body fluids (semen), and cell suspensions containing $\leq 1.0 \times 10^6$ cells as well as lysates derived from Proteinase K digested cells and tissues.

1. Add 4 volumes of **Genomic Lysis Buffer** to each volume of liquid sample. (For example, for 200 μ l of sample, add 800 μ l of **Genomic Lysis Buffer**). Mix briefly by vortexing 4-6 seconds.

Note: For Proteinase K digested material, add 4 volumes of **Genomic Lysis Buffer** to each volume of lysate then mix briefly by vortexing. Centrifuge the mixture at top speed ($>10,000 \times g$) for 5 minutes. Transfer up to 1.0 ml supernatant to a new microcentrifuge tube and proceed to Step 2.

2. Ensure the **ZymoBeads™** are fully resuspended by vortexing. Add 10 μ l **ZymoBeads™** to the supernatant. Mix by inversion, and then incubate at room temperature for 5 minutes. Centrifuge the tube at $1,500 \times g$ for 1 minute. Discard the supernatant.
3. Add 200 μ l of **Genomic Lysis Buffer** to the **ZymoBeads™**. Resuspend pellet by pipetting up and down. Centrifuge at $1,500 \times g$ for 1 minute. Discard the supernatant.
4. Add 200 μ l of **DNA Pre-Wash Buffer** to the **ZymoBeads™**. Resuspend the pellet and then centrifuge at $1,500 \times g$ for 1 minute. Discard the supernatant.
5. Add 500 μ l of **g-DNA Wash Buffer** to the **ZymoBeads™**. Resuspend the pellet and then centrifuge at $1,500 \times g$ for 1 minute. Discard the supernatant. Recentrifuge briefly and remove any residual wash buffer.
6. Add $\geq 35 \mu$ l of **Elution Buffer** or water, resuspend pellet by pipetting up and down, and then centrifuge at $10,000 \times g$ for 1 minute.
7. Collect the supernatant. The supernatant contains purified DNA that can be used immediately or stored (e.g., -20°C) for later use.

Troubleshooting:

1. **DNA degradation:** Check for DNase contamination. All reagents and components supplied with the **ZR Genomic DNA I Kit™** are DNase-free. However, DNase contamination could result during the processing of some samples. Check pipets, pipet tips, microcentrifuge tubes, etc., and exercise the appropriate precautions during the DNA purification procedure.
2. **DNA is not performing well in subsequent experiments:** Ensure the correct volume of **Genomic Lysis Buffer** has been added to the sample. Also, make sure all centrifugation steps are completed for the indicated times and speeds (rcfs). Failure to do so may result in incomplete washing, which may cause salts to be eluted with the DNA affecting quantitation and subsequent experiments including enzymatic processes like PCR.
3. **RNA contamination:** The buffers and spin columns provided in this kit are designed to efficiently remove RNA during the DNA purification procedure. However, additional RNA removal (e.g., digestion with RNase A) may be necessary for subsequent applications sensitive to trace amounts of RNA.

Ordering Information

Product Description	Catalog No.	Kit Size
ZR Genomic DNA I Kit™	D3004	100 purifications
ZR Genomic DNA I Kit™	D3005	400 purifications

For Individual Sale	Catalog No.	Amount
Genomic Lysis Buffer	D3004-1-50	50 ml
	D3004-1-100	100 ml
DNA Pre-Wash Buffer	D3004-5-15	15 ml
	D3004-5-30	30 ml
	D3004-5-50	50 ml
g-DNA Wash Buffer	D3004-2-50	50 ml
	D3004-2-100	100 ml
DNA Elution Buffer	D3004-4-4	4 ml
	D3004-4-16	16 ml
ZymoBeads™	D3004-3-1	1 ml
	D3004-3-4	4 x 1ml

Popular DNA Purification & Analysis Products from Zymo Research

Product	Description	Kit Size (Preps)	Catalog No. (column format)
DNA Clean & Concentrator-5™	Clean & concentrate DNA from any reaction or “crude” preparation in 2 minutes. A 6 µl minimum elution volume allows for highly concentrated DNA. Designed for samples containing up to 5 µg of DNA.	50	D4003 (uncapped)
		200	D4004 (uncapped)
		50	D4013 (capped)
		200	D4014 (capped)
DNA Clean & Concentrator-25™	Clean & concentrate DNA in minutes. 25 µl minimum elution volume allows for highly concentrated DNA. Designed for purifying up to 25 µg of DNA.	50	D4005 (uncapped)
		200	D4006 (uncapped)
		50	D4033 (capped)
		200	D4034 (capped)
DNA Clean & Concentrator-100™	Clean & concentrate DNA in minutes. 100 µl minimum elution volume allows for highly concentrated DNA. Designed for purifying up to 100 µg of DNA.	25	D4029
		50	D4030
DNA Clean & Concentrator-500™	Clean & concentrate DNA in minutes. 1 ml minimum elution volume allows for highly concentrated DNA. Designed for samples containing up to 500 µg of DNA.	10	D4031
		20	D4032
ZR-96 DNA Clean & Concentrator-5™	Quick (15 minute), high-output recovery of pure DNA from PCR, endonuclease digestions, plasmid preparations, etc. 10-15 µl minimum elution volume allows for highly concentrated DNA. Designed for samples containing up to 5 µg of DNA.	2x96	D4023
		4x96	D4024
Zymoclean™ Gel DNA Recovery Kit	Purify DNA from high and low-melting agarose gels in minutes	50	D4001
		200	D4002
ZR-96 Zymoclean™ Gel DNA Recovery Kit	High-throughput DNA purification from high and low-melting agarose gels.	2x96	D4021
		4x96	D4022
Pinpoint Slide DNA Isolation System™	Recover genomic DNA from paraffin-embedded or fresh tissue sections for PCR. Ideal for isolating DNA from clinical tissue samples.	50	D3001
Zyppy™ Plasmid Miniprep Kit	Pellet-Free™ plasmid DNA purification in minutes: (alkaline lysis/spin column format for low 30 µl elution volume).	50	D4036
		100	D4019
		400	D4020
Zyppy™ Plasmid Midiprep Kit	Pellet-Free™ plasmid DNA purification in minutes: (alkaline lysis/spin column format and 150 µl minimum elution volume).	25	D4025
		50	D4026
Zyppy™ Plasmid Maxiprep Kit	High-purity plasmid DNA purification in minutes: (alkaline lysis/spin column format and 2 ml minimum elution volume).	10	D4027
		20	D4028
ZR Genomic DNA I Kit™	Genomic DNA isolation from whole blood, tissue culture cells, solid tissue and liquid samples. (Silica bead format is scalable to fit your requirements).	100	D3004
		400	D3005
ZR Genomic DNA II Kit™	Genomic DNA purification from whole blood, tissue culture cells, solid tissue and liquid samples. No requirement for beads or phenol chloroform.	50	D3006 (uncapped)
		200	D3007 (uncapped)
		50	D3024 (capped)
		200	D3025 (capped)
ZR-96 Genomic DNA Kit™	High-output genomic DNA purification from whole blood, tissue culture cells, solid tissue and liquid samples. No requirement for beads or phenol chloroform.	2x96	D3010
		4x96	D3011
ZR Soil Microbe DNA Kit™	Simple, rapid isolation of humic-free, PCR-quality genomic DNA from soil microbes.	50	D6001
ZR Fungal/Bacterial DNA Kit™	Simple, rapid isolation of PCR-quality genomic DNA from fungi.	50	D6005
ZR Fecal DNA Kit™	Simple, rapid isolation of PCR-quality genomic DNA from feces.	50	D6010
ZR Viral DNA Kit™	Isolation of viral DNA from cell-free body fluids or sample mixtures containing cells at concentrations less than 10 ⁵ cells per ml.	50	D3015
		200	D3016
ZR-96 Viral DNA Kit™	High-output (96-well) isolation of viral DNA from cell-free body fluids or sample mixtures containing cells at concentrations less than 10 ⁵ cells per ml.	2x96	D3017
		4x96	D3018
EZ DNA Methylation™ Kit	Streamlined kit for the conversion of unmethylated cytosines in DNA to uracil via the <u>chemical-denaturation</u> of DNA using our specially designed CT Conversion Reagent™. DNA is then desulphonated and subsequently cleaned using <i>Fast-Spin</i> column technology. Ultra-pure recovered DNA can be used for PCR and bisulfite sequencing applications.	50	D5001
		200	D5002
		2x96	D5003 (Shallow-well)
		2x96	D5004 (Deep-well)
EZ DNA Methylation-Gold™ Kit	Streamlined kit for the conversion of unmethylated cytosines in DNA to uracil via <u>heat-denaturation</u> of DNA using our specially designed CT Conversion Reagent™. DNA is then desulphonated and subsequently cleaned using <i>Fast-Spin</i> column technology. Ultra-pure recovered DNA can be used for PCR and bisulfite sequencing applications. <i>3 hour processing time!</i>	50	D5005
		200	D5006
		2x96	D5007 (Shallow-well)
		2x96	D5008 (Deep-well)

*Bulk quantities are available upon request. Please contact: busdev@zymoresearch.com or call 1-888-882-9682 for assistance.

ZYMO RESEARCH CORP.

Toll Free: 1-888-882-9682 • Fax: 1-714-288-9643 • Web: www.zymoresearch.com • E-mail: info@zymoresearch.com