



**ZYMO RESEARCH**

*The Beauty of Science is to Make Things Simple*

# INSTRUCTION MANUAL

## ZR Genomic DNA II Kit™

Catalog Nos. **D3006, D3007, D3024, & D3025**

### Highlights

- Easy purification of high quality DNA from whole blood, plasma, serum, body fluids, buffy coat, lymphocytes, tissue, swabs or cultured cells in less than 15 minutes using innovative *Fast-Spin* column technology.
- Compatible with commonly used anticoagulants (i.e., EDTA, heparin, citrate).
- Excludes the use of Proteinase K and organic denaturants.
- Isolated DNA is ideal for PCR, endonuclease digestion, bisulfite conversion/methylation detection, sequencing, genotyping, etc.

### Contents

Product Contents & Specifications.....	1
Product Description .....	2
Protocols	
Whole Blood, Serum, and Plasma Samples..	3
Buccal Cell Samples.....	4
Solid Tissue Samples.....	4
Cell Monolayer Samples.....	5
Cell Suspensions and Proteinase K Digested Samples.....	6
Troubleshooting.....	7
Ordering Information.....	7
List of DNA-Related Products.....	8

Satisfaction of all Zymo Research products is guaranteed. If you are not satisfied with this product please call 1-888-882-9682.

## Product Contents

ZR Genomic DNA II Kit™ (Kit Size)	D3006 (uncapped) D3024 (capped) (50 Preps.)	D3007 (uncapped) D3025 (capped) (200 Preps.)	Storage Temperature
<b>Genomic Lysis Buffer</b>	50 ml	2 x 100 ml	Room Temp.
<b>DNA Pre-Wash Buffer</b>	15 ml	50 ml	Room Temp.
<b>g-DNA Wash Buffer</b>	50 ml	100 ml	Room Temp.
<b>DNA Elution Buffer</b>	4 ml	10 ml	Room Temp.
<b>Zymo-Spin™ Columns*</b>	50 columns	200 columns	Room Temp.
<b>Collection Tubes</b>	100 tubes	400 tubes	Room Temp.
<b>Instruction Manual</b>	1	1	-

Note - Integrity of kit components is guaranteed for up to one year from date of purchase. Reagents are routinely tested on a lot-to-lot basis to ensure they provide maximal performance and reliability.

\* Kits D3006 and D3007 are supplied with uncapped columns (Zymo-Spin™ IIN), whereas kits D3024 and D3025 are supplied with capped columns (Zymo-Spin™ IIC).

## Specifications

- **Sample Sources** – Whole blood, plasma, or serum from humans, mice, rats, etc. Also, tissue, cells from culture, buccal cells, as well as a variety of biological liquids are effectively processed using this kit.
- **DNA Purity** – High-quality DNA is eluted with **DNA Elution Buffer** or water. DNA is especially well suited for PCR and other downstream applications. Typical absorption indices are  $A_{260}/A_{280} > 1.8$ .
- **DNA Size Limits** – Capable of recovering genomic DNA up to and above 40 kb. In most instances, mitochondrial DNA and viral DNA (if present) will also be recovered.
- **DNA Recovery** – Typically, up to 25 µg total DNA is eluted into  $\geq 50$  µl (30 µl minimum) **DNA Elution Buffer** or water. Human whole blood will yield 3 - 7 µg DNA per 100 µl blood sampled. Mammalian tissues yield: 1 - 3 µg DNA per mg skeletal, heart, and brain tissues and 3 - 5 µg DNA per mg liver, kidney and lung tissues.
- **Product Detergent Tolerance** –  $\leq 5\%$  Triton X-100,  $\leq 5\%$  Tween-20,  $\leq 5\%$  Sarkosyl,  $\leq 0.1\%$  SDS.
- **Equipment** – microcentrifuge, vortex


Note - ™ Trademarks of Zymo Research Corporation. This product is for research use only and should only be used by trained professionals. Some reagents included with this kit are irritants. Wear protective gloves and eye protection. Follow the safety guidelines and rules enacted by your research institution or facility.

ZYMO RESEARCH CORP.

Toll Free: 1-888-882-9682 • Fax: 1-714-288-9643 • Web: [www.zymoresearch.com](http://www.zymoresearch.com) • E-mail: [info@zymoresearch.com](mailto:info@zymoresearch.com)

## Product Description

The **ZR Genomic DNA II Kit™** is a simple procedure for the rapid isolation of total DNA (e.g., genomic, mitochondrial, viral) from a variety of biological sample sources. This product has been optimized for maximal recovery of ultra-pure DNA without RNA contamination and is compatible with whole blood (fresh or stored), serum, plasma, buffy coat, solid tissue, bone marrow and buccal cells, cells from culture, and many biological liquid samples. For processing, simply add the specially formulated **Genomic Lysis Buffer** to a sample, vortex, and transfer the mixture to the supplied **Zymo-Spin™ Column**. There is no need for organic denaturants or Proteinase K digestion because of the unique chemistries featured in the kit. Instead, the product features *Fast-Spin* column technology to yield high-quality, purified DNA in just minutes (see below). PCR inhibitors are effectively removed during the purification process. DNA purified using the **ZR Genomic DNA II Kit™** is suitable for PCR, nucleotide blotting, DNA sequencing, endonuclease digestion, bisulfite conversion/methylation analysis, and other downstream applications.



High yield/quality DNA is successfully isolated from porcine whole blood using the **ZR Genomic DNA II Kit™**. Equivalent amounts (100 µl) of blood were processed without Proteinase K using the **ZR Genomic DNA II Kit™** in half the time as compared to the kits from suppliers Q and S. Equal volumes of eluted DNA were then analyzed (in duplicate) in a 0.8% (w/v) TAE/agarose/ethidium bromide gel. The size marker "M" is a 1 kb ladder (Zymo Research).

The **ZR Genomic DNA I Kit™ (D3004, D3005)** and **ZR Serum DNA Kit™ (D3013)** are recommended for yields >25 µg per prep. Both feature silica beads instead of a spin column.

Zymo Research offers the **EZ DNA Methylation™ (D5001, D5002, D5003)**, **EZ DNA Methylation-Gold™ (D5005, D5006, D5007, D5008)** and **EZ DNA Methylation-Direct™ (D5020, D5021, D5022, D5023)** Kits for rapid, precise DNA methylation detection.

For **Technical Assistance**, please contact 1-888-882-9682 or E-mail tech@zymoresearch.com.

## **PROTOCOLS**

### **WHOLE BLOOD, SERUM, AND PLASMA SAMPLES**

The following is for the purification of DNA from 50 µl whole blood, serum or plasma (the volumes can be adjusted depending on your requirements). Fresh, frozen, or preserved blood (in EDTA, citrate, or heparin) can be used. If material cannot be processed immediately, the sample can be “stabilized” for later processing (as noted below) although the immediate processing of blood samples is recommended.

1. Add 200 µl of **Genomic Lysis Buffer** to 50 µl of blood, serum, or plasma. Mix completely by vortexing 4 - 6 seconds, then let stand 5 minutes at room temperature.

**Note:** Add 200 µl Genomic Lysis Buffer to all samples < 50 µl. For samples larger than 50 µl, add a proportional amount of Genomic Lysis Buffer (e.g., for 100 µl blood, add 400 µl Genomic Lysis Buffer).

2. Transfer the mixture to a **Zymo-Spin™ Column** in a **Collection Tube**. Centrifuge at 10,000 x g for one minute. Discard the **Collection Tube** with the flow through.
3. Transfer the **Zymo-Spin™ Column** to a new **Collection Tube**. Add 200 µl of **DNA Pre-Wash Buffer** to the spin column. Centrifuge at 10,000 x g for one minute.
4. Add 500 µl of **g-DNA Wash Buffer** to the spin column. Centrifuge at 10,000 x g for one minute.
5. Transfer the spin column to a clean microcentrifuge tube. Add ≥ 50 µl **DNA Elution Buffer** or water to the spin column. Incubate 2 - 5 minutes at room temperature, then centrifuge at top speed for 30 seconds to elute the DNA. The eluted DNA can be used immediately for molecular based applications or stored ≤ -20°C for future use.

**Delayed Processing of Blood Samples:** The immediate processing of blood with this kit is recommended. However, if blood cannot be processed immediately, samples can be “stabilized” in **Genomic Lysis Buffer** for processing at a later time. To do this, add four volumes of **Genomic Lysis Buffer** to each volume of whole blood (4:1), then vortex. Blood samples mixed with **Genomic Lysis Buffer** can be stored at room temperature for 7 days, 0 - 4°C for up to 15 days, -20°C for 2 months, or < -70°C for many years. Samples stored at ≤ 4°C should reach room temperature prior to processing. Begin at Step 2 in the standard protocol (above) when purifying DNA from blood samples stabilized in **Genomic Lysis Buffer**.

The column capacity is ~1 ml.

Elution of DNA from the column is dependent on pH and temperature. If water is used, ensure the pH is > 6.0. Also, the total yield may be improved by eluting the DNA with Elution Buffer or water pre-equilibrated to 60-70°C.

## **BUCCAL CELL SAMPLES**

Buccal cells can be isolated using a rinse- or swab-based isolation method.

- A. **Rinse Method:** Vigorously rinse 10 - 20 ml of saline solution or mouthwash orally for 30 seconds. The more vigorous the rinsing action, the more cells that will be recovered. Spit the saline into a 50 ml tube and pellet the cells at 1,500 rpm for 5 minutes. Discard the supernatant without disturbing the cell pellet. Add 500  $\mu$ l of **Genomic Lysis Buffer** to the pellet then vortex 4 - 6 seconds, then let stand at room temperature for 5 minutes.
- B. **Swab Isolation Method:** Thoroughly rinse mouth out before isolating cells. Brush the inside of the cheek with a buccal swab for 15 seconds (approximately 20 brushes), making sure to cover the entire area of the inner cheek. Rinse the brush into a microcentrifuge tube using 500  $\mu$ l of **Genomic Lysis Buffer** then vortex 4 - 6 seconds, then let stand at room temperature for 5 minutes.
  1. Transfer the mixture to a **Zymo-Spin™ Column** in a **Collection Tube**. Centrifuge at 10,000 x g for one minute. Discard the **Collection Tube** with the flow through.
  2. Transfer the **Zymo-Spin™ Column** to a new **Collection Tube**. Add 200  $\mu$ l of **DNA Pre-Wash Buffer** to the spin column. Centrifuge at 10,000 x g for one minute.
  3. Add 500  $\mu$ l of **g-DNA Wash Buffer** to the spin column. Centrifuge at 10,000 x g for one minute.
  4. Transfer the spin column to a clean microcentrifuge tube. Add  $\geq$  50  $\mu$ l **DNA Elution Buffer** or water to the spin column. Incubate 2 - 5 minutes at room temperature, then centrifuge at top speed for 30 seconds to elute the DNA. The eluted DNA can be used immediately for molecular based applications or stored  $\leq$  -20°C for future use.

## **SOLID TISSUE SAMPLES**

**Note:** For Proteinase K digested materials (e.g., tailsnips) follow the protocol for **CELL SUSPENSIONS AND PROTEINASE K DIGESTED SAMPLES** (pg. 6). Otherwise...mechanically homogenize up to 20 mg of fresh or frozen tissue in 500  $\mu$ l of **Genomic Lysis Buffer**. Increase the reagents proportionally if more than 20 mg of solid tissue is used or reduce the reagents to half the volume described in the protocol should the tissue sampled be less than 5 mg.

1. Centrifuge the lysate at top speed (10,000 x g) for 5 minutes. Making sure not to disturb the pelleted debris, transfer the supernatant to a **Zymo-Spin™ Column** in a **Collection Tube** and centrifuge at 10,000 x g for one minute. Discard the **Collection Tube** with the flow through.
2. Transfer the **Zymo-Spin™ Column** to a new **Collection Tube**. Add 200  $\mu$ l of **DNA Pre-Wash Buffer** to the spin column. Centrifuge at 10,000 x g for one minute.

Soft tissue samples are readily homogenized using our **Squisher™-Single**, **Squisher™-8**, and **Squisher™-96** products.

Typical yields are: 1 - 3  $\mu$ g DNA per mg skeletal, heart, and brain tissues and 3 - 5  $\mu$ g per mg liver, kidney, and lung tissues.

3. Add 500  $\mu\text{l}$  of **g-DNA Wash Buffer** to the spin column. Centrifuge at 10,000  $\times g$  for one minute.
4. Transfer the spin column to a clean microcentrifuge tube. Add  $\geq 50 \mu\text{l}$  **DNA Elution Buffer** or water to the spin column. Incubate 2 - 5 minutes at room temperature, then centrifuge at top speed for 30 seconds to elute the DNA. The eluted DNA can be used immediately for molecular based applications or stored  $\leq -20 \text{ }^\circ\text{C}$  for future use.

### **CELL MONOLAYER SAMPLES**

Generally, no more than  $5 \times 10^6$  cells should be sampled, for larger samples will exceed the binding capacity of the spin column.

It may be necessary to centrifuge the sample mixture before transferring the supernatant to the **Zymo-Spin™ Column** to remove particulate matter that may clog the column.

The column capacity is  $\sim 1 \text{ ml}$ .

The following procedure is designed for up to  $1.0 \times 10^6$  monolayer cells (roughly equal to one well of a 6-well plate or  $\frac{1}{2}$  of a T25 flask). Although cell types and culture conditions may vary, the protocol will work with high-density growth cells (e.g., HeLa cells) as well as with low-density growth cells (e.g., neuronal cells). The procedure may be scaled up or down for increases or decreases in the amounts of monolayer cells sampled (see the **Guidelines for Monolayer Cell DNA Isolation** below).

1. Trypsinize or manually scrape adherent cells from the growth surface of a culture flask or plate. Centrifuge the cell suspension at approximately 500  $\times g$  for 5 minutes. Remove the supernatant and add 500  $\mu\text{l}$  of **Genomic Lysis Buffer** directly to the pellet. Resuspend pellet by vortexing 4 - 6 seconds and let stand for 5 minutes at room temperature.
2. Transfer the mixture to a **Zymo-Spin™ Column** in a **Collection Tube**. Centrifuge at 10,000  $\times g$  for one minute. Discard the **Collection Tube** with the flow through.
3. Transfer the **Zymo-Spin™ Column** to a new **Collection Tube**. Add 200  $\mu\text{l}$  of **DNA Pre-Wash Buffer** to the spin column. Centrifuge at 10,000  $\times g$  for one minute.
4. Add 500  $\mu\text{l}$  of **g-DNA Wash Buffer** to the spin column. Centrifuge at 10,000  $\times g$  for one minute.
5. Transfer the spin column to a clean microcentrifuge tube. Add  $\geq 50 \mu\text{l}$  **DNA Elution Buffer** or water to the spin column. Incubate 2 - 5 minutes at room temperature, then centrifuge at top speed for 30 seconds to elute the DNA. The eluted DNA can be used immediately for molecular based applications or stored  $\leq -20 \text{ }^\circ\text{C}$  for future use.

**Guidelines for Monolayer Cell DNA Isolation:** The above procedure is designed for the processing of  $0.1 - 1.0 \times 10^6$  cells. However, cell numbers (growth densities) can vary between different cell types. Table 1 (pg. 6) provides an approximation of what can be recovered from different culture containers for high-density growth cells like CV1 and HeLa cells. The table provides a guideline when adjusting reagent volumes used in the protocol. If the cell number is less than  $1.0 \times 10^5$ , use half the volume of **Genomic Lysis Buffer** given in the standard protocol (above). If the cell number is greater than  $1.0 \times 10^6$ , use double the volume of **Genomic Lysis Buffer**.



**Table 1: Culture Plate/Flask Growth Area (cm<sup>2</sup>) and Cell Number**

Culture Container	Well /Flask Surface Area	Cell Number
96-well plate	0.32 - 0.6 cm <sup>2</sup>	4 - 5 x 10 <sup>4</sup>
24-well plate	2 cm <sup>2</sup>	1 - 3 x 10 <sup>5</sup>
12-well plate	4 cm <sup>2</sup>	4 - 5 x 10 <sup>5</sup>
6 -well plate	9.5 cm <sup>2</sup>	0.5 - 1 x 10 <sup>6</sup>
T25 Culture Flask	25 cm <sup>2</sup>	2 - 3 x 10 <sup>6</sup>
T75 Culture Flask	75 cm <sup>2</sup>	0.6 - 1 x 10 <sup>7</sup>
T175 Culture Flask	175 cm <sup>2</sup>	2 - 3 x 10 <sup>7</sup>

### **CELL SUSPENSIONS AND PROTEINASE K DIGESTED SAMPLES**

The following protocol is designed for up to 200 µl of biological liquid sample including CSF, buffy coat, body fluids (semen), and cell suspensions containing less than 5.0 x 10<sup>6</sup> cells as well as lysates derived from Proteinase K digested samples.

1. Add 4 volumes of **Genomic Lysis Buffer** to each volume of liquid sample. (For example, for 200 µl of sample, add 800 µl of **Genomic Lysis Buffer**). Mix briefly by vortexing, then let stand at room temperature for 5 minutes.

**Note:** For Proteinase K digested material, add 4 volumes of **Genomic Lysis Buffer** to each volume of lysate then mix briefly by vortexing. Centrifuge the mixture at top speed (>10,000 x g) for 5 minutes. Transfer up to 1.0 ml supernatant to the Zymo-Spin™ Column in Step 2.

2. Transfer the mixture to a **Zymo-Spin™ Column** in a **Collection Tube**. Centrifuge at 10,000 x g for one minute. Discard the **Collection Tube** with the flow through.
3. Transfer the **Zymo-Spin™ Column** to a new **Collection Tube**. Add 200 µl of **DNA Pre-Wash Buffer** to the spin column. Centrifuge at 10,000 x g for one minute.
4. Add 500 µl of **g-DNA Wash Buffer** to the spin column. Centrifuge at 10,000 x g for one minute.
5. Transfer the spin column to a clean microcentrifuge tube. Add ≥ 50 µl **DNA Elution Buffer** or water to the spin column. Incubate 2 - 5 minutes at room temperature, then centrifuge at top speed for 30 seconds to elute the DNA. The eluted DNA can be used immediately for molecular based applications or stored ≤ -20°C for future use.

Cells should be processed directly from biological fluids or from suspension in PBS, TE, or compatible buffers.

The column capacity is ~1 ml.

Typical yields from Proteinase K digested tissues are: 1-3 µg DNA per mg skeletal, heart, and brain tissues and 3-5 µg per mg liver, kidney, and lung tissues.

**Troubleshooting:**

1. **DNA degradation:** Check for DNase contamination. All reagents and components supplied with the **ZR Genomic DNA II Kit™** are DNase-free. However, DNase contamination could result during the processing of some samples. Check pipets, pipet tips, microcentrifuge tubes, etc., and exercise the appropriate precautions during the DNA purification procedure.
2. **DNA is not performing well in subsequent experiments:** Ensure the correct volume of **Genomic Lysis Buffer** has been added to the sample. Also, make sure all centrifugation steps are completed for the indicated times and speeds (rcfs). Failure to do so may result in incomplete washing, which may cause salts to be eluted with the DNA affecting quantitation and subsequent experiments including enzymatic processes like PCR.
3. **RNA contamination:** The buffers and spin columns provided in this kit are designed to efficiently remove RNA during the DNA purification procedure. However, additional RNA removal (e.g., digestion with RNase A) may be necessary for subsequent applications sensitive to trace amounts of RNA.

**Ordering Information**

Product Description	Catalog No.	Kit Size
<b>ZR Genomic DNA II Kit™</b> supplied w/ <b>uncapped</b> columns	D3006	50 preps.
<b>ZR Genomic DNA II Kit™</b> supplied w/ <b>uncapped</b> columns	D3007	200 preps.
<b>ZR Genomic DNA II Kit™</b> supplied w/ <b>capped</b> columns	D3024	50 preps.
<b>ZR Genomic DNA II Kit™</b> supplied w/ <b>capped</b> columns	D3025	200 preps.

For Individual Sale	Catalog No.	Amount
<b>Genomic Lysis Buffer</b>	D3004-1-50	50 ml
	D3004-1-100	100 ml
<b>DNA Pre-Wash Buffer</b>	D3004-5-15	15 ml
	D3004-5-30	30 ml
	D3004-5-50	50 ml
<b>g-DNA Wash Buffer</b>	D3004-2-50	50 ml
	D3004-2-100	100 ml
<b>DNA Elution Buffer</b>	D3004-4-4	4 ml
	D3004-4-10	10 ml
<b>Zymo-Spin™ IIN Columns (uncapped)</b>	C1019-50	50 columns
	C1019-250	250 columns
<b>Zymo-Spin™ IIC Columns (capped)</b>	C1011-50	50 columns
	C1011-250	250 columns
<b>Collection Tubes</b>	C1001-50	50 tubes
	C1001-500	500 tubes
	C1001-1000	1,000 tubes

ZYMO RESEARCH CORP.

Toll Free: 1-888-882-9682 • Fax: 1-714-288-9643 • Web: [www.zymoresearch.com](http://www.zymoresearch.com) • E-mail: [info@zymoresearch.com](mailto:info@zymoresearch.com)



## Popular DNA Purification & Analysis Products from Zymo Research

Product	Description	Kit Size (Preps)	Catalog No. (column format)
<b>DNA Clean &amp; Concentrator-5™</b>	Clean & concentrate DNA from any reaction or “crude” preparation in 2 minutes. A 6 µl minimum elution volume allows for highly concentrated DNA. Designed for samples containing up to 5 µg of DNA.	50	<b>D4003</b> (uncapped)
		200	<b>D4004</b> (uncapped)
		50	<b>D4013</b> (capped)
		200	<b>D4014</b> (capped)
<b>DNA Clean &amp; Concentrator-25™</b>	Clean & concentrate DNA in minutes. 25 µl minimum elution volume allows for highly concentrated DNA. Designed for purifying up to 25 µg of DNA.	50	<b>D4005</b> (uncapped)
		200	<b>D4006</b> (uncapped)
		50	<b>D4033</b> (capped)
		200	<b>D4034</b> (capped)
<b>DNA Clean &amp; Concentrator-100™</b>	Clean & concentrate DNA in minutes. 100 µl minimum elution volume allows for highly concentrated DNA. Designed for purifying up to 100 µg of DNA.	25	<b>D4029</b>
		50	<b>D4030</b>
<b>DNA Clean &amp; Concentrator-500™</b>	Clean & concentrate DNA in minutes. 1 ml minimum elution volume allows for highly concentrated DNA. Designed for samples containing up to 500 µg of DNA.	10	<b>D4031</b>
		20	<b>D4032</b>
<b>ZR-96 DNA Clean &amp; Concentrator-5™</b>	Quick (15 minute), high-output recovery of pure DNA from PCR, endonuclease digestions, plasmid preparations, etc. 10-15 µl minimum elution volume allows for highly concentrated DNA. Designed for samples containing up to 5 µg of DNA.	2x96	<b>D4023</b>
		4x96	<b>D4024</b>
<b>Zymoclean™ Gel DNA Recovery Kit</b>	Purify DNA from high and low-melting agarose gels in minutes	50 200	<b>D4001</b> <b>D4002</b>
<b>ZR-96 Zymoclean™ Gel DNA Recovery Kit</b>	High-throughput DNA purification from high and low-melting agarose gels.	2x96	<b>D4021</b>
		4x96	<b>D4022</b>
<b>Pinpoint Slide DNA Isolation System™</b>	Recover genomic DNA from paraffin-embedded or fresh tissue sections for PCR. Ideal for isolating DNA from clinical tissue samples.	50	<b>D3001</b>
<b>Zyppy™ Plasmid Miniprep Kit</b>	Pellet-Free™ plasmid DNA purification in minutes: (alkaline lysis/spin column format for low 30 µl elution volume).	50	<b>D4036</b>
		100	<b>D4019</b>
		400	<b>D4020</b>
<b>Zyppy™ Plasmid Midiprep Kit</b>	Pellet-Free™ plasmid DNA purification in minutes: (alkaline lysis/spin column format and 150 µl minimum elution volume).	25	<b>D4025</b>
		50	<b>D4026</b>
<b>Zyppy™ Plasmid Maxiprep Kit</b>	High-purity plasmid DNA purification in minutes: (alkaline lysis/spin column format and 2 ml minimum elution volume).	10	<b>D4027</b>
		20	<b>D4028</b>
<b>ZR Genomic DNA I Kit™</b>	Genomic DNA isolation from whole blood, tissue culture cells, solid tissue and liquid samples. (Silica bead format is scalable to fit your requirements).	100	<b>D3004</b>
		400	<b>D3005</b>
<b>ZR Genomic DNA II Kit™</b>	Genomic DNA purification from whole blood, tissue culture cells, solid tissue and liquid samples. No requirement for beads or phenol chloroform.	50	<b>D3006</b> (uncapped)
		200	<b>D3007</b> (uncapped)
		50	<b>D3024</b> (capped)
		200	<b>D3025</b> (capped)
<b>ZR-96 Genomic DNA Kit™</b>	High-output genomic DNA purification from whole blood, tissue culture cells, solid tissue and liquid samples. No requirement for beads or phenol chloroform.	2x96	<b>D3010</b>
		4x96	<b>D3011</b>
<b>ZR Soil Microbe DNA Kit™</b>	Simple, rapid isolation of humic-free, PCR-quality genomic DNA from soil microbes.	50	<b>D6001</b>
<b>ZR Fungal/Bacterial DNA Kit™</b>	Simple, rapid isolation of PCR-quality genomic DNA from fungi.	50	<b>D6005</b>
<b>ZR Fecal DNA Kit™</b>	Simple, rapid isolation of PCR-quality genomic DNA from feces.	50	<b>D6010</b>
<b>ZR Viral DNA Kit™</b>	Isolation of viral DNA from cell-free body fluids or sample mixtures containing cells at concentrations less than 10 <sup>5</sup> cells per ml.	50	<b>D3015</b>
		200	<b>D3016</b>
<b>ZR-96 Viral DNA Kit™</b>	High-output (96-well) isolation of viral DNA from cell-free body fluids or sample mixtures containing cells at concentrations less than 10 <sup>5</sup> cells per ml.	2x96	<b>D3017</b>
		4x96	<b>D3018</b>
<b>EZ DNA Methylation™ Kit</b>	Streamlined kit for the conversion of unmethylated cytosines in DNA to uracil via the <u>chemical-denaturation</u> of DNA using our specially designed CT Conversion Reagent™. DNA is then desulphonated and subsequently cleaned using <i>Fast-Spin</i> column technology. Ultra-pure recovered DNA can be used for PCR and bisulfite sequencing applications.	50	<b>D5001</b>
		200	<b>D5002</b>
		2x96	<b>D5003</b> (Shallow-well)
		2x96	<b>D5004</b> (Deep-well)
<b>EZ DNA Methylation-Gold™ Kit</b>	Streamlined kit for the conversion of unmethylated cytosines in DNA to uracil via <u>heat-denaturation</u> of DNA using our specially designed CT Conversion Reagent™. DNA is then desulphonated and subsequently cleaned using <i>Fast-Spin</i> column technology. Ultra-pure recovered DNA can be used for PCR and bisulfite sequencing applications. <i>3 hour processing time!</i>	50	<b>D5005</b>
		200	<b>D5006</b>
		2x96	<b>D5007</b> (Shallow-well)
		2x96	<b>D5008</b> (Deep-well)

\*Bulk quantities are available upon request. Please contact: [busdev@zymoresearch.com](mailto:busdev@zymoresearch.com) or call 1-888-882-9682 for assistance.

ZYMO RESEARCH CORP.

Toll Free: 1-888-882-9682 • Fax: 1-714-288-9643 • Web: [www.zymoresearch.com](http://www.zymoresearch.com) • E-mail: [info@zymoresearch.com](mailto:info@zymoresearch.com)