



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

INSTRUCTION MANUAL

ZR-96 RNA Clean & Concentrator™

Catalog No. R1080

Highlights

- Quick (*20 minute*) high-throughput (96-well) method for cleaning and concentrating RNA.
- *Fast-Spin* plate technology allows RNA to be eluted into minimal volumes ($\geq 10 \mu\text{l}$).
- Eluted RNA is ultra clean and ready for subsequent analysis and molecular manipulation.

Contents

Product Contents.....	1
Specifications.....	1
Product Description.....	2
Buffer Preparation.....	3
Protocol.....	3
Ordering Information.....	4
Appendix.....	5
Related Products.....	6

Satisfaction of all Zymo Research products is guaranteed. If you should be dissatisfied with this product please call 1-888-882-9682.

Product Contents

ZR-96 RNA Clean & Concentrator™ (Kit Size)	R1080 (2x 96 Preps.)	Storage Temperature
RNA Binding Buffer	100 ml	Room Temp.
RNA Prep Buffer	2x 25 ml	Room Temp.
RNA Wash Buffer¹ (concentrate)	24 ml	Room Temp.
DNase/RNase-Free Water	4 ml	Room Temp.
Zymo-Spin™ I-96 Plate	2	Room Temp.
Collection Plate	2	Room Temp.
Elution Plate	2	Room Temp.
Cover Foil	4	Room Temp.
Instruction Manual	1	-

Note - Integrity of kit components is guaranteed for up to one year from date of purchase. Reagents are routinely tested on a lot-to-lot basis to ensure they provide the highest performance and reliability.

¹ Add 96 ml 100% ethanol (104 ml of 95% ethanol) to the 24 ml **RNA Wash Buffer** concentrate before use.

Specifications

- **Sample Sources** – RNA fragments (≥ 17 nucleotides); e.g., DNase treated RNA, *in vitro* transcription products.
- **Format** – 96-well Plate
- **RNA Purity** – High quality RNA ($A_{260}/A_{280} > 1.8$, $A_{260}/A_{230} > 1.8$) suitable for all downstream RNA-based manipulations.
- **RNA Recovery** – Typically, RNA is eluted into as little as ≥ 10 μ l RNase-free water allowing for a highly concentrated sample. The RNA binding capacity of the supplied **Zymo-Spin™ I-96 Plate** is ~ 25 μ g/well.
- **RNA Storage** – RNA is eluted with RNase-free water and can be stored at ≤ -70 °C. The addition of RNase inhibitors is optional but highly recommended for prolonged storage.
- **Equipment Needed** – Microcentrifuge

Note:

For DNA clean-up see the **ZR-96 DNA Clean & Concentrator™-5** (Catalog Nos. D4023, D4024).

Note - ™ Trademarks of Zymo Research Corporation. This product is for research use only and should only be used by trained professionals. Some reagents included with this kit are irritants. Wear protective gloves and eye protection. Follow the safety guidelines and rules enacted by your research institution or facility.

ZYMO RESEARCH CORP.

Toll Free: 1-888-882-9682 • Fax: 1-714-288-9643 • Web: www.zymoresearch.com • E-mail: info@zymoresearch.com

Product Description


The **ZR-96 RNA Clean & Concentrator™** provides a simple and reliable method for large-scale (96-well) purification and concentration of up to ~25 µg (per well) of high-quality RT-PCR-ready RNA. The procedure is based on the use of a unique single-buffer system and *Fast-Spin* plate technology.

The procedure features a single-buffer system that allows for efficient RNA adsorption from up to 96 samples onto the matrices of the supplied **Zymo-Spin™ I-96 Plate**. The RNA is washed twice then eluted and concentrated into ≥10 µl of RNase-free water.

The kit allows for small RNA recovery or removal by alternating the first two steps of the purification procedure. RNA fragments (≥17 bases) can be safely treated and recovered using this kit. The result is highly-concentrated, purified RNA that is suitable for subsequent RNA-based methods including RT-PCR, hybridization, etc.

The entire RNA isolation procedure typically takes about 20 minutes.

For **Assistance**, please contact Zymo Research Technical Support at 1-888-882-9682 or e-mail tech@zymoresearch.com.



- ✓ Durable Construction
- ✓ Deep-Well Design
- ✓ Small Elution Volume

The **Zymo-Spin™ I-96 Plate** design.

Make sure guidelines are followed to ensure the RNA isolation procedure is performed in an RNase-free environment.

Notes:

¹ RNA species ≥ 17 nt (e.g., miRNAs) will be recovered.

² Add 100 μ l **RNA Binding Buffer** (Step 1a) or 0.5 \times **RNA Binding Buffer** (Step 1b) to all samples ≤ 50 μ l.

³ **Example:** Mix by pipetting or vortex a sealed 96-well plate. **Collection Plate** may be also used for the sample preparation. (If a **Collection Plate** is used for mixing, use the **Cover Foil** and quick spin prior to the next step.)

⁴ The sample capacity of each well of the **Collection Plate** is ~ 800 μ l. It may be necessary to load and spin the **Zymo-Spin™ I-96 Plate** multiple times to process samples > 800 μ l.

⁵ Small RNA species ≤ 200 nt will be removed

⁶ **Optional DNase treatment:** Following Step 3, samples can be in-column DNase treated. See **Appendix** (page 5).

⁷ For quantitative recovery of RNA samples, increase the elution volume (> 15 μ l) and/or repeat the elution Step 7.

⁸ To prevent an evaporation and to store purified RNA samples on the **Elution Plate** use the **Cover Foil**.

Buffer Preparation

Before starting, add 96 ml of 100% ethanol (104 ml of 95% ethanol) to the 24 ml **RNA Wash Buffer** concentrate.

Protocol

Follow part A or B of the protocol before proceeding to Step 3.

A. To recover total RNA including small RNAs¹:

- 1a. Mix 1 volume of an RNA sample with 2 volumes of the **RNA Binding Buffer** (e.g., 100 μ l RNA sample and 200 μ l buffer)² by first aliquoting the **RNA Binding Buffer** to the wells of the **Zymo-Spin™ I-96 Plate** mounted on a **Collection Plate** and then adding RNA samples. Mix well^{3,4}.
- 2a. Add 1 volume 95-100% ethanol to the mixture from Step 1a (e.g., 300 μ l ethanol and 300 μ l mixture), mix well and proceed to Step 3 below.

B. To remove the small RNA species from an RNA sample⁵:

- 1b. Prepare a 0.5 \times **RNA Binding Buffer** by mixing 1 volume of **RNA Binding Buffer** with 1 volume ethanol (95-100%).
- 2b. Mix 1 volume of an RNA sample with 2 volumes of the 0.5 \times **RNA Binding Buffer** from Step 1b (e.g., 100 μ l RNA sample and 200 μ l buffer)² by first aliquoting the 0.5 \times **RNA Binding Buffer** to the wells of the **Zymo-Spin™ I-96 Plate** mounted on a **Collection Plate** and then adding RNA samples^{3,4}.

3. Centrifuge the **Zymo-Spin™ I-96 Plate** from Step 2a or 2b, mounted on a **Collection Plate**, at 5,000 $\times g$ for 5 minutes⁶. Discard the flow-through and mount the **Zymo-Spin™ I-96 Plate** back onto the **Collection Plate**.
4. Add 400 μ l **RNA Prep Buffer** to each well of the **Zymo-Spin™ I-96 Plate**. Centrifuge at 5,000 $\times g$ for 5 minutes. Discard the flow-through.
5. Add 800 μ l **RNA Wash Buffer** to each well of the **Zymo-Spin™ I-96 Plate**. Centrifuge at 5,000 $\times g$ for 2 minutes. Discard the flow-through. Repeat the wash step with 400 μ l **RNA Wash Buffer** and centrifuge at 5,000 $\times g$ for 2 minutes. Discard the flow-through.
6. Centrifuge the **Zymo-Spin™ I-96 Plate** at 5,000 $\times g$ for 5 minutes in the emptied **Collection Plate** to ensure complete removal of the wash buffer.
7. Transfer the **Zymo-Spin™ I-96 Plate** onto the **Elution Plate**. Add ≥ 10 μ l **RNase-Free Water** directly to the column matrix in each well⁷ and let stand for 1 minute at room temperature. To elute the RNA, centrifuge at 5,000 $\times g$ for 5 minutes. RNA can be used immediately or stored⁸ at ≤ -70 °C.

Ordering Information

Product Description	Catalog No.	Kit Size
ZR-96 RNA Clean & Concentrator™	R1080	2x 96 Preps.

For Individual Sale	Catalog No.	Amount
RNA Binding Buffer	R1013-2-50	50 ml
	R1013-2-100	100 ml
RNA Prep Buffer	R1060-2-10	10 ml
	R1060-2-25	25 ml
RNA Wash Buffer (concentrate)	R1003-3-6	6 ml
	R1003-3-12	12 ml
	R1003-3-24	24 ml
	R1003-3-48	48 ml
DNase/RNase-Free Water	W1001-1	1 ml
	W1001-4	4 ml
	W1001-6	6 ml
	W1001-10	10 ml
Zymo-Spin™ I-96 Plates	C2004	2 plates
Collection Plates	C2002	2 plates
Elution Plates	C2003	2 plates

Fast-Spin plate technology employed in the ZR RNA isolation products efficiently removes the majority of DNA during RNA purification and is satisfactory for most RNA-based applications. However, if necessary, complete removal of DNA can be achieved by performing a DNase I digestion.

Appendix

In-Column DNase Digestion

The DNase digestion procedure can be performed using any source of RNase-free DNase I together with its 10x reaction buffer (e.g., 100 U **RNase-free DNase I (1 U/μl) w/ 10x Reaction Buffer** – Zymo Research Cat. No. **E1007**). DNase I maintains activity in the **RNA Wash Buffer** provided in this kit.

1. Make the following DNase I cocktail (for each sample/well to be treated):

RNase-Free DNase I	3 μl (1 U/μl)
10x Reaction Buffer	3 μl
RNA Wash Buffer	24 μl

2. Following Step 3 of the RNA isolation protocol¹, add 400 μl **RNA Wash Buffer** to the **Zymo-Spin™ I-96 Plate** mounted on a **Collection Plate** and centrifuge at 5,000 x *g* for 5 minutes. Discard the flow through.
3. Add 30 μl DNase I cocktail from Step 1 above directly to the matrix of the **Zymo-Spin™ I-96 Plate**. Keep the **Zymo-Spin™ I-96 Plate** mounted on a **Collection Plate**.
4. Incubate the plate at 25-37°C for ≥15 minutes², then centrifuge at 5,000 x *g* for 5 minutes.
5. Continue with Step 4 of the RNA isolation protocol³.

Notes:

¹ See page 3, Protocol – step 3.

² The temperature optimum for DNase I activity is at 37 °C. An optimal incubation time may vary.

³ See page 3, Protocol – step 4.

Related Products

Product	Description	Kit Size	Cat No.
<i>RNA Clean-up, Concentration & Recovery</i>			
RNA Clean & Concentrator™-5	Spin Column Format (up to 5 µg/prep.)	50 Preps. 200 Preps.	R1015 R1016
RNA Clean & Concentrator™-25	Spin Column Format (up to 25 µg/prep.)	50 Preps. 200 Preps.	R1017 R1018
ZR-96 RNA Clean & Concentrator™	96-Well Format (up to 25 µg/well)	2x96 Preps.	R1080
DNA-Free RNA Kit™	Spin Column Format (up to 5 µg/prep.)	50 Preps. 200 Preps.	R1013 R1014
Zymoclean™ Gel RNA Recovery Kit	Spin Column Format (up to 5 µg/prep.)	50 Preps.	R1011
ZR small-RNA™ PAGE Recovery Kit	Spin Column Format (up to 5 µg/prep.)	20 Preps.	R1070