

General Protocol

Immunofluorescence Kit

For Research Use Only

Peninsula Laboratories, Inc.
A Member of the Bachem Group
305 Old County Rd., San Carlos, CA 94070
Tel +1 800 922 1516
+1 650 592 5392
Fax +1 650 595 4071

Table of Contents

Introduction.....	2
Materials Provided	2
Materials Required but not Provided	2
Preparation Procedure	3
Staining Procedure	3
Summary of Immunofluorescence Technique	4
Technical Notes.....	4
References.....	5
Ordering Information	6
Technical Support	6
Guarantee and Limitation of Remedy	6

Introduction

This kit is designed to stain a specific peptide antigen in fixed tissue by an indirect immunofluorescence technique. The indirect immunofluorescence technique uses the primary antibody directed against the specific antigen (peptide) in the tissue followed by a secondary antibody, which binds to the primary antibody and is conjugated with a fluorescent dye, such as fluorescein-isothiocyanate (FITC). Positive staining will show an emission of fluorescent light (FITC: green to yellow) in a dark field. We recommend the use of frozen tissue sections. The kit provides sufficient reagents for 50 slides.

Materials Provided

- Primary antibody: peptide-specific rabbit (or guinea pig) antiserum, 10 ml (1:200) (lyophilized powder).
- Secondary antibody: fluorescent dye-conjugated goat antibody directed against rabbit (or guinea pig) immunoglobulins, 100 µl.
- Negative control: normal rabbit (or guinea pig) serum, 5 µl (lyophilized).
- Phosphate buffered saline (PBS) concentrate (5x), 10 ml, pH 7.4.
- 3% Triton X-100/water, 5 ml.
- Mounting media: Glycerin-PBS (9:1), 2 ml.
- Normal goat serum (NGS, 1:10), 13 ml (lyophilized powder).

Materials Required but not Provided

- Rinsing buffer: 0.1 M phosphate buffered saline (PBS), pH 7.4, containing: 19 mM sodium phosphate monobasic, monohydrate; 81 mM sodium phosphate dibasic, anhydrous; 0.05 M sodium chloride. Store at 4 °C.
- Glass microscope slides coated with gelatin; cover slips.
- Incubation chamber or dishes.
- Cryostat.
- O.C.T. (Tissue-Tek, Miles Laboratories, Inc.) embedding medium.
- Dark field microscope (using transmitted or reflected light) with filters which are appropriate for FITC detection (excitation at 440-490 nm, emission at 520-560 nm).

Preparation Procedures

- Rehydrate 0.1 M PBS with 45 ml of distilled water. Add 5 ml of 3% Triton X-100/water to the PBS solution. This will give 50 ml of 0.1 M PBS containing 0.3% Triton X-100.
- Rehydrate primary antiserum with 10 ml of 0.1 M PBS containing 0.3% Triton X-100. This will give 10 ml of primary antiserum (1:200), which is sufficient for 50 slides using 200 µl of antiserum per slide.
- Rehydrate NGS with 13 ml of 0.1 M PBS containing 0.3% Triton X-100. This will give 13 ml of NGS (1:10).
- Rehydrate normal rabbit serum with 1 ml of 0.1M PBS. This will give 1 ml of NRS at a concentration of 1:200.
- Dilute the fluorescent secondary antibody with 9.9 ml of 0.1 M PBS containing 0.3% Triton X-100. This will give 10 ml of the fluorescent secondary antibody (1:100).

Staining Procedure

Day 1


- Bring the frozen tissue sections to room temperature and defrost for 5 minutes.
- Wash slides 3 times for 10 minutes each with 0.1 M PBS.
- Blot excess PBS from the slides with an absorbent wipe (**Do not touch the tissue or allow it to dry**).
- Apply 200 µl of NGS (1:10) on the tissue sections. Incubate the slides at room temperature for 30 minutes in a closed incubation chamber.
- Blot slides (**Do not touch the tissue or allow it to dry**).
- Apply approximately 200 µl of the primary antiserum (1:200) or negative control on tissue sections. Incubate the slides at 4 °C for 18-24 hours in the same chamber.

Day 2

- Wash the slides 3 times for 10 minutes each with 0.1 M PBS.
- Blot the excess PBS from the slides with an absorbent wipe (**Do not touch the tissue or allow it to dry**).
- Apply 200 µl of the fluorescent secondary antibody (either goat anti-rabbit or anti-guinea pig IgG) (1:100). Incubate the slides at room temperature for 60 minutes in the same chamber.
- Wash the slides 3 times for 10 minutes each with 0.1 M PBS.
- Drain and blot the excess PBS (**Do not touch the tissue or allow it to dry**). Coverslip with glycerin-PBS (9:1).

- F) Examine the stained tissue sections under a dark field microscope for the presence of a specific fluorescent staining pattern.

Summary of Immunofluorescence Technique



This procedure is not intended to replace the preparation steps, precautions and procedures shown in this protocol. Please consult those sections for complete directions.

Technical Notes

- A) Storage: 4 °C.
B) Expiration date: 3 months.
C) Tissue fixation: 4% paraformaldehyde in 0.1 M PBS, pH 7.4.
D) Dilution: The 1:200 dilution of the primary antibody is suitable for most antisera. However, if the background is too high or if the staining is either too strong or unclear, the antibody can be diluted further.
E) Crossreactivity: Use the same procedure as in the negative control. The concentrations of the tested peptides are 1-10 µM.

References

1. A.H. Coons
Fluorescent antibody methods. General cytochemical methods. p. 399-422, ed. J. Danielli. Academic Press, New York (1958)
2. M. Schultzberg et al.
Distribution of peptide- and catecholamine-containing neurons in the gastro-intestinal tract of rat and guinea-pig: immunohistochemical studies with antisera to substance P, vasoactive intestinal polypeptide, enkephalins, somatostatin, gastrin/cholecystokinin, neurotensin and dopamine beta-hydroxylase. *Neuroscience* **5**, 689-744 (1980)
3. K. Jessen et al.
Immunohistochemical studies of the enteric nervous system in tissue culture and in situ: localization of vasoactive intestinal polypeptide (VIP), substance-P and enkephalin immunoreactive nerves in the guinea-pig gut. *Neuroscience* **5**, 1717-1735 (1980)
4. J. Haaijman
Labeling of proteins with fluorescent dyes: quantitative aspects of immunofluorescence microscopy. *Immunohistochemistry*. p. 47-85, ed. A. Cuello. John Wiley & Sons (1984)
5. Y.N. Wang et al.
Immunocytochemical localization of neuropeptide Y-like immunoreactivity in adrenergic and non-adrenergic neurons of the rat gastrointestinal tract. *Peptides* **8**, 145-151 (1987)
6. R.M. Kobayashi et al.
Regional distribution of neurotensin and somatostatin in rat brain. *Brain Res.* **126**, 584-588 (1977)
7. R.N. Burgus et al.
Primary structure of somatostatin, a hypothalamic peptide that inhibits the secretion of pituitary growth hormone. *Proc. Natl. Acad. Sci. USA* **70**, 684-688 (1973)

Ordering Information

For additional kits or our most current catalog products, please visit our web site at www.bachem.com or call us at:

Tel +1 800 922 1516

+1 650 592 5392

Fax +1 650 595 4071

Technical Support

If you need technical information or assistance with assay procedures, please call our Technical Service Department at (800) 922-1516 or (650) 592-5392. Our staff will be happy to answer your questions about this or any other products.

Guarantee and Limitation of Remedy

The Peninsula Laboratories, Inc. makes no guarantee of any kind, expressed or implied, which extends beyond the description of the material in this kit, except that these materials and this kit will meet our specifications at the time of delivery. Customer's remedy and Peninsula Laboratories, Inc.'s sole liabilities hereunder are limited at Peninsula Laboratories, Inc.'s option to either refund the purchase price or to replace the material that does not meet our specifications.

By the acceptance of our products, the customer indemnifies and holds Peninsula Laboratories, Inc. harmless against, and assumes all liability for, the consequences of its use or misuse by the customer, its employees or others.