TERESA

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
Fig.1. Serum hIFN α concentration at various time points after intramuscular gene delivery and electroporation

乙肝治疗性 DNA 疫苗恒河猴药效试验 (Yong-Gang Zhao, Baowei Peng, Yuhong Xu, *et al.* Anti-HBV immune responses in rhesus macaques elicited by electroporation mediated DNA vaccination. *Vaccine*, 2006, 24:897-903)


Fig.2. HBV S protein specific antibody titers in rhesus macaques during and after three treatments of DNA vaccination with or without electroporation (n = 3, mean \pm S.E.M.). Peripheral blood samples were drawn at specific times and sera were isolated and assayed for antibody titers using the ELISA Kit; The three treatment times were indicated using arrows pointing at the X-axis. (■) 0.1 mg dose with electroporation; (●) 0.5 mg dose with electroporation; (▲) 1 mg dose with electroporation; (▼) 1 mg dose without electroporation. The dotted line in the bottom of the figure indicates the threshold of 10 mIU/mL of anti-HBVs antibody

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1. F.-Q. Yang, Y.-Y. Yu, G.-Q. Wang, *et al.* A pilot randomized controlled trial of dual-plasmid HBV DNA vaccine mediated by *in vivo* electroporation in chronic hepatitis B patients under lamivudine chemotherapy. *J Viral Hepatitis*, 2012, DOI: 10.1111/j.1365-2893.2012.01589.x
2. Teng Zuo, Xuanling Shi, Linqi Zhang, *et al.* Comprehensive analysis of pathogen-specific antibody response *in vivo* based on an antigen library displayed on surface of yeast. *J Biol Chem*, 2011, 286(38): 33511-33519.
3. Yong-Gang Zhao, Yuhong Xu. Electroporation-mediated HBV DNA vaccination in primate models. *Methods Mol Biol*, 2008, 423(6): 487-495.
4. Jian-hui NIE, Lin-qi ZHANG, You-chun WANG, *et al.* Comparison of the Immunogenicities of HIV-1 mutants based on structural modification of *env*. *Virologica Sinica*, 2008, 23(4): 233-246.
5. Baowei Peng, Yonggang Zhao, Yuhong Xu, *et al.* Electric pulses applied prior to intramuscular DNA vaccination greatly improve the vaccine immunogenicity. *Vaccine*, 2007, 25: 2064 - 73.
6. Yong-gang ZHAO, Hui-li LU, Yu-hong XU, *et al.* Inhibitory effect of Ca²⁺ on *in vivo* gene transfer by electroporation. *Acta Pharmacologica Sinica*, 2006, 27 (3): 307 - 310.
7. Yong-Gang Zhao, Baowei Peng, Yuhong Xu, *et al.* Anti-HBV immune responses in rhesus macaques elicited by electroporation mediated DNA vaccination. *Vaccine*, 2006, 24: 897 - 903.
8. Pang WK, Zhao Y, Xu Y, *et al.* Therapeutic protein production *in vivo* after electroporation-assisted intramuscular gene delivery. *Conf Proc IEEE Eng Med Biol Soc*, 2006, 6: 6229-32.
9. Baowei Peng, Yonggang Zhao, Yuhong Xu, *et al.* *In vivo* plasmid DNA electroporation resulted in transfection of satellite cells and lasting transgene expression in regenerated muscle fibers. *Biochem Biophys Res Commun*, 2005, 338(3):1490 - 8.

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Shanghai Teresa Healthcare Sci-Tech Co., Ltd.