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	Direct UV A280 Absorbance	Colorimetric	Fluorescence
Sample purity	Samples must be purified as contaminants may interfere with measurement.	No purification necessary.	No purification necessary.
Buffer compatibility	Buffers with strong UV absorbance may be unsuitable (see next panel).	Some assays are sensitive to detergents or reducing reagents, which can artificially perturb or enhance color development.	Some assays are sensitive to buffers with primary amines or detergents, which perturb fluorescent signal.
Other considerations	Knowledge of an E1% value or molecular weight and molar extinction coefficient are required to calculate mg/mL concentration.	Colorimetric signals vary between proteins, therefore standards must be carefully chosen in order to minimize differences in signal between the standard and sample proteins.	Fluorescent assays typically have a lower detection limit than colorimetric assays or direct UV A280 measurement, but are limited by the maximum measurable protein concentration (fig. 2).

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	Direct UV A280 Absorbance	Colorimetric	Fluorescence
Preparation	None	Production of working solutions is required in some assays; fresh standards need to be made for a calibration curve.	Production of working solutions is required in most assays; fresh standards need to be made for a calibration curve.
Incubation	None	Time required to stabilize colorimetric signal is between 10 and 60 minutes, depending on the assay used (fig. 5).	Typical time required to stabilize fluorescent signal is 10 – 30 minutes, depending on the assay used.
Measurement	~ 5 seconds using a NanoDrop	Colorimetric assays require more time than the direct UV A280 method as a calibration curve must be established prior to unknown sample quantification.	Fluorescent assays require more time than the direct UV A280 method as a calibration curve must be established prior to unknown sample quantification.

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©2012 赛默飞世尔科技公司保留所有权利。所有商标均归赛默飞世尔科技公司及其子公司所有。