

BIOSCIENCES[™] G-Biosciences, St Louis, MO, USA ♦ 1-800-628-7730 ♦ 1-314-991-6034 ♦ technical@GBiosciences.com

BLOT-QuickBlocker TM

INTRODUCTION

Protein solutions are used for blocking excess binding sites in western blots and other immunochemical applications. Each blocking solution will impart certain characteristics on the assay system. BLOT-QuickBlocker[™] is a modified milk protein, and is good for blocking western blot membranes.

ITEM(S) SUPPLIED	Cat # 786-011
<i>BLOT</i> -QuickBlocker [™]	175g

STORAGE CONDITIONS

It is shipped at ambient temperature. Store at RT, stable for 1 year.

PROTOCOL

- Prepare 5% solution of *BLOT*-QuickBlocker[™] in TBS /TBST or PBS/PBST. 1.
- Add 15-20ml of 5% solution of *BLOT*-QuickBlocker^{$^{\text{TM}}$} on to the blotting membrane. Make sure the membrane is 2. submerged in the solution.
- 3. Incubate the membrane with BLOT-QuickBlockerTM for 30-60 minutes or more at room temperature with gentle
- Proceed with the blot detection protocol as per your specific application. 4.

RELATED PRODUCTS

- 1. Protease Arrest (Cat # 786-108) A cocktail of protease inhibitors for use during protein extraction and purification. Protease ArrestTM inhibits a broad spectrum of serine, cysteine and metalloprotease as well as calpains.
- 2. NI-Protein Assay (Cat # 786-005) A protein assay that is free from interference of common laboratory agents, including reducing agents, detergents, dyes, EDTA, etc.
- 3. Foldase (Cat # 786-185) A kit for optimization of protein folding procedure. Includes a selection of Foldase, i.e. agents to facilitate protein folding.
- 4. Strategy (Cat # 786-184) A kit for optimization of purification of novel proteins. Includes a selection of agents, tools and simple to follow protocols.

<u>NOTE</u>: For other related products, visit our web site at <u>www.GBiosciences.com</u> or contact us.

