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Website: www.bioer.com.cn

at 12,000 $\times g$. Discard the flow-through.

9. Repeat step 8th once.
10. Centrifuge for an additional 1 minute at 12,000 $\times g$ and transfer the Spin column to a sterile 1.5ml micro centrifuge tube. Recommend to centrifuge according to this step; otherwise, there will be residual liquid in the column
11. Add 50 μ l Elution Buffer, ddH₂O or TE Buffer to the Spin column and let it stand for 1 minute at room temperature.
The volume of elution buffer could be adjusted according to needs.
12. Centrifuge for 1 minute at 12,000 $\times g$. The buffer in the micro centrifuge tube contains the plasmid DNA.
13. The purified plasmid DNA can be used directly for kinds of downstream molecular biological experiments. Store at -20 $^{\circ}$ C if not used immediately.

Troubleshooting

Get on plasmid DNA

If the plasmid DNA is not found in elution buffer, please check whether the ethanol had been added to Wash Buffer according to the volume be marked on bottle label

Low plasmid DNA yields

- 1) Please make sure that the bacteria was cultured in the right way,
- 2) The bacteria cells should be resuspended completely.
- 3) Incubate the Elution Buffer in 30~60 $^{\circ}$ C, it will increase the yields.

Absorbance problem

- 1) Absorbance is the difference from sample and criterion, please use the Elution Buffer to adjust to zero value and dilute the sample.
- 2) If the ratio of OD₂₆₀/ OD₂₃₀ is low, wash the spin column for one more time.
- 3) In the case of low ratio of OD₂₆₀₋₃₂₀/ OD₂₈₀₋₃₂₀, there is protein contamination, please add Neutralization Buffer, and then centrifuge buffer with sufficient rotating speed, thus to make precipitation compact; be careful to pipette supernatant so as to avoid pipette precipitation.
- 4) If the ratio of OD₂₆₀₋₃₂₀/ OD₂₈₀₋₃₂₀ is high, add Rnase A to Resuspension Buffer (100 μ g/ml)

Electrophoresis problem

- 1) If there is genomic DNA in the result. Please invert the tube gently(step4 and 5).
- 2) If there is RNA in the result, add Rnase A to Resuspension Buffer (100 μ g/ml).

DNA Analysis

⊗ Absorbance analysis

Get some plasmid DNA, diluted in a advisable factor with elution buffer.

Survey the OD₂₆₀, OD₂₈₀ and OD₃₂₀

expressions: concentration (μ g/ml) = 50 \times OD₂₆₀ \times dilution fact


target: $2.0 \geq OD_{260-320} / OD_{280-320} \geq 1.8$

Notice: $1.0 \geq OD_{260} \geq 0.1$, the result of ratio is much reliable.

⊗ Agarose Gel Analysis

0.8~1% Agarose gel


Example 1: plasmid DNA electrophoresis

Example 2: enzymatic reactions analysis

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Important notes

1. The RNase solution should be all added into the Resuspension Buffer before use, mix and store at 2-8°C.
2. Add ethanol (as the volume be marked on bottle label) to Wash Buffer and mix well.
3. If the Lysis Buffer and Neutralization Buffer precipitated, it should be Redissolve by warming to 37°C. please not vortex Lysis Buffer acutely.
4. Please close the lid immediately after using Lysis Buffer so as to avoid acidification.
5. The Kit can extract high-quality plasmid DNA from 1-5ml E.coli overnight cultured.
6. The suitable volume is 50ul for Elution Buffer, user can adjust its volume if necessary.

Procedure

1. **Add 1-1.5ml cultured bacteria to 1.5ml micro centrifuge tube.**
2. **centrifuge at 10, 000rpm (8,000~10,000×g) for 30s, and discard the supernatant.**
Step 1th and 2th could be repeated for more than one time to collect enough cells.
3. **Resuspend pelleted bacterial cells in 250µl Resuspension Buffer and No cell clumps should be visible after resuspension of the pellets.**
4. **Add 250µl Lysis Buffer and gently invert the tube 4-6 times to mix.**
Do not vortex, as this will result in shearing of genomic DNA. Do not allow this step for more than 5 minutes.
5. **Add 350µl Neutralization Buffer and gently invert the tube 4-6 times to mix.**
The solution should become cloudy and no local precipitate should be visible.
6. **Centrifuge for 10minutes at 13,000rpm (>14,000×g) until a compact white pellet form.**
7. **Apply the supernatant to the Spin column and centrifuge for 30-60 seconds at 6,000×g. Discard the flow-through.**
8. **Add 650µl Wash Buffer to the Spin column and Centrifuge for 30-60 seconds**

Kit Components(100T)

Component	Amount
Resuspension Buffer	25ml
Lysis Buffer	25ml
Neutralization Buffer	35ml
Wash Buffer	60ml
Elution Buffer	20ml
RNase solution	1
Spin columns	100
Handbook	1copy

Storage and Transportation

- ◆ The Kit should be stored at room temperature(15~25°C), but the RNase solution should be stored at 2~8°C. The kit can be stored for up to 18 months by this method. After addition of RNase solution, Resuspension Buffer should be stored 2~8°C.
- ◆ The kit can be transported at room temperature.

Introduction

The Kit provides a fast, simple, and cost-effective plasmid miniprep method for routine molecular biology laboratory applications. Plasmid DNA can be purified from 1–5ml of overnight cultures of *E. coli*. The DNA isolated by this Kit is ready for downstream applications such as restriction enzyme digestion, sequencing, PCR/Real-time PCR and other downstream experiments.

Apparatus and Materials to Be Supplied by the User

- * sterile 1.5ml micro centrifuge tubes
- * 10µl/100µl/1000µl tips
- * microcentrifuge capable of 14,000 × g
- * Absolute ethanol
- * Vortex mixer

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