



# Cancer MicroRNA Array

Catalog Number AP-0003

(For Research Use Only)

## Introduction

Newly discovered microRNAs (miRNAs) are important to the regulation of gene expression (1). Up to 30% of mammalian genes might be regulated by miRNAs. So far, more than 400 miRNAs have been identified in the human genome and many of them differ only by one or a few nucleotides. Expression of mature miRNAs is tissue-specific (2) and the abundance of miRNAs varies in several orders of magnitude (3). More importantly, misregulation of miRNA expression might contribute to human cancers (4). Systematic profiling of miRNA expression displays unique signatures in a number of cancers (5).

Based on its proprietary technology, Signosis has developed a highly sensitive and discriminative array technology for monitoring miRNA. It combines oligo-ligation assay-based detection and T7 transcription-based linear amplification, providing a highly sensitive and specific array assay. Via ligation, selection, and transcription of a pair of oligos that are hybridized to each specific miRNA, the targets are amplified and the array monitored. By assigning unique tag sequences to individual isoforms, they can be easily differentiated. The whole procedure is simple and straightforward. We are currently offering an array that targets 132 of the most well studied miRNAs, each with cited literature, which is able to facilitate the comparison and discovery of differentiated miRNA expression in different samples.

## Principles of the technology

miRNAs are different from large messenger RNAs in three aspects; (1) miRNAs are small sized molecules with quite a big difference in abundance, (2) mature miRNAs co-exist with their precursor pre-miRNA and pri-miRNA, only differing in length, and (3) many miRNAs are very closely related in sequences, such as isoforms, differing by only one or a few nucleotides. Therefore, the conventional microarray technologies cannot directly be applied to analyzing these molecules. A number of miRNA microarray products are commercially available, but they are either tedious in requiring pre-isolation of the microRNA, lack discriminative power to differentiate between isoforms, or are not sensitive enough to monitor low abundant miRNAs.

In our array assay, each miRNA molecule is targeted by two oligos, each that hybridizes a half molecule of the

target miRNA to form a RNA/DNA duplex. When the sequences are perfectly matched, they are aligned with the miRNA and the joint can be ligated by DNA ligase (figure 1). A single nucleotide difference among miRNAs will block either the hybridization or the ligation, so that miRNA isoform can be differentiated. Due to the small size of miRNA, the hybrid might not be stable; therefore we introduce the stacking sequences. By extending these two oligos along with their complementary oligos the stability is increased. Once the pair of oligos are ligated, the ligated molecules are subjected to linear amplification via T7 transcription into RNA in the presence of biotin-UTP, which are used as probes for array hybridization. To differentiate each isoform, we assigned unique tag sequences to the ligation oligos, so that single nucleotide differences are converted into unique tag sequences. Therefore, each isoform can be easily distinguished by array hybridization.

We offer the miRNA profiling assay kit to profile the expression of 132 cancer-associated miRNAs and their isoforms. The procedure is simple and straight forward, including three steps: (1) mix total RNA with provided oligos to form miRNA/oligo hybrids; (2) select the hybrids and remove free oligos, and ligate miRNA-directed pairing of oligos to become a single DNA; and (3) Amplify the ligated DNA with T7 transcription.




Figure 1: Diagram of miRNA array assay

## Materials provided with the kit

30 µl Cancer oligo Mix (-20°C)  
8 µl Array Detection Oligo (-20°C)  
400 µl Annealing buffer (RT)  
15 µl Magnetic streptavidin beads (4 °C)  
120 µl Beads binding buffer (RT)  
1 ml Bead wash buffer (RT)  
6 µl Ligase (-20°C)  
250 µl Ligation buffer (-20°C)  
60 µl Extension mix (-20°C)  
60 µl Labeling mix (-20°C)  
6 µl T7 RNA polymerase (-20°C)  
30ml 1x Hybridization buffer (RT)  
30ml 5x Hybridization wash buffer (RT)  
60ml Blocking buffer (RT)  
50 µl Streptavidin-HRP conjugate (4 °C)  
40ml 5x Detection wash buffer (RT)  
3.6 ml Substrate A (4 °C)  
3.6 ml Substrate B (4 °C)  
3 Array membranes (RT)  
3 Detection sheets (RT)

## Material required but not provided

Magnetic stand (96 well plate)  
RNase-free water  
PCR machine  
Hybridization oven  
Washing tray  
50ml Centrifuge tubes (Corning tubes are recommended, cat#430290) or hybridization bottles  
0.2ml PCR tube  
Alpha Innotech image or equivalent image system or X-ray film

## Sample and reagent preparation before starting experiment

- Dilute the 5x Hybridization wash buffer and 5x Detection wash buffer to 1x buffer:  
1x Hybridization wash buffer:  
30ml 5x Hybridization washing buffer  
120ml ddH<sub>2</sub>O  
1x Detection wash buffer:  
40ml 5x Detection washing buffer  
160 ml ddH<sub>2</sub>O
- Prewarm 1x Hybridization buffer, 1x Hybridization wash buffer at 42°C for 1 hour or until the buffers are clear without visible precipitation before using.
- Pre-hybridization can be done during T7 RNA transcription at the Step 4.

## Assay Procedure

### 1. Annealing of miRNA with Oligo mix

- (1) Sample preparation
- X µl 5ug total RNA or 10ng isolated miRNA
  - 10 µl Cancer oligo mix
  - 2 µl Array Detection Oligo
  - 20 µl Annealing buffer
  - X µl ddH<sub>2</sub>O

-----  
40ul

- (2) Incubate on the PCR machine at 72°C for 5 minutes and 53°C for 90 minutes.

### 2. Selection of miRNA/oligo hybrids

- (1) Beads washing
- Resuspend the beads by gently tapping the tube to obtain a homogeneous suspension.
  - Transfer 5 µl beads to a 0.2ml PCR tube (the size of the tube should fit into the magnetic stand.
  - Add 100 µl of annealing buffer to the tube and then place onto the magnetic stand for 30 seconds.
  - Aspirate out the liquid.
  - Remove the tube from the magnetic stand.

- (2) Beads selection
- Add 40 µl of Bead binding buffer to 40 µl annealed miRNA/oligo hybrid mix from Step 1, transfer to the washed beads from Step 2(1) and resuspend the beads in the solution.
  - Incubate at 37°C for 30 minutes.
  - Place the bead mixture on the magnetic stand for 30 seconds, and aspirate out the buffer. The beads will remain on the side of tube.
  - Remove the tube from the magnetic stand and add 100 µl of Bead wash buffer to resuspend the beads, pipette gently up and down, and put the tube on the magnetic stand for 30 seconds, aspirate the buffer.
  - Repeat the washing step once.

### 3. Ligation of miRNA-directed oligos to form a single molecule

- (1) Add 50 µl of Ligation buffer to resuspend the beads, pipette gently up and down, place the tube on the magnetic stand for 30 seconds, and aspirate the buffer.
- (2) Remove the tube from the magnetic stand. Add 20 µl of ligation buffer to resuspend the beads, then add 2µl ligase to the resuspended beads and incubate at 37° for 90 minutes.

#### 4. T7 RNA transcription of ligated molecule

- (1) Add 100  $\mu$ l of Bead washing buffer directly to the 20  $\mu$ l ligation reaction mix from Step 3, place the tube on the magnetic stand for 30 seconds, and aspirate the buffer.
- (2) Remove the tube from the magnetic stand and add 20  $\mu$ l of Extension mix to resuspend the beads.
- (3) Incubate the mixture on the PCR machine at 94°C for 2 minutes, 54°C for 1 minute, 72°C for 1.5 minutes, and 94°C for 30 seconds.
- (4) Place the reaction tube on the magnetic stand for 30 seconds. **Immediately** transfer the 20  $\mu$ l of the extension mix to a fresh tube (keep the solution and toss the beads).
- (5) Add 20  $\mu$ l of Labeling mix and 1  $\mu$ l of T7 RNA polymerase to the tube.
- (6) Incubate the mixture at 37°C for 1 hour.
- (7) The transcribed RNA is ready for hybridization.

#### 5. Pre-hybridization and hybridization


- (1) Place each array membrane into a 50 ml tube. Wet the membrane by filling the tube with dH<sub>2</sub>O, then carefully decant the water. The side of the membrane with the spotted oligos should face into the middle of the tube.
- (2) Add 4 ml of prewarmed 1x Hybridization buffer to each tube. Incubate the tubes in a hybridization oven at 42°C for at least 30-60 minutes.
- (3) Decant 1x Hybridization buffer, replace with 4ml of prewarmed 1x Hybridization buffer and 40  $\mu$ l of transcribed RNA to the prehybridized membrane and incubate overnight in a hybridization oven at 42°C.
- (4) Decant the hybridization mixture from each bottle and wash each membrane as follows:
  - Rinse the membrane with 20 ml Hybridization washing buffer, and decant liquid.
  - Incubate the membrane with 20 ml Hybridization wash buffer at 42°C for 20 minutes. Decant liquid.

#### 6. Detection

- (1) Using forceps, carefully transfer the membrane from the hybridization tube to a container (an empty 200  $\mu$ l pipette tips box). Each box could have two membranes, one at each side of the box.
- (2) Rinse the membrane with 10 ml of 1X Detection wash buffer.
- (3) Block the membrane with 15 ml of Blocking buffer for 30 minutes at room temperature with moderate shaking.
- (4) Dilute 15  $\mu$ l of Streptavidin-HRP conjugate with 1 ml of the 1X Blocking buffer and transfer to the container. **Do not** add the HRP diluted solution directly onto the membrane.
- (5) Continue shaking the membrane for 45 min at room temperature.
- (6) Decant the Blocking buffer and wash three times at room temperature with 15 ml of 1x Detection washing buffer, 10 minutes each wash.

- (7) Mix equal amounts of Substrate A and B. Place the membrane on the bottom side of detection sheet on a flat surface and overlay the membrane with 2 ml of substrate solution. To ensure that the solution remains evenly distributed over the membrane when enveloped by the detection sheet: gently lower the top side of the detection sheet halfway over the membrane then pull back up slightly to allow the solution to flow back over the membrane. Then slowly lay the top sheet down completely without trapping air bubbles. Incubate at room temperature for 5 minutes.
- (8) Remove excess substrate by gently applying pressure over the top sheet using a paper towel. Expose the membranes using either Hyperfilm ECL (2-10 min) or a chemiluminescence imaging system (i.e., FluorChem imager from Alpha Innotech). With either method, experiment with different exposure times.
- (9) Use the schematic diagram of human miRNA array I to identify the spots on the array.

#### Example of cancer miRNA array analysis.



5 $\mu$ g of HeGP2, 293 and HeLa total RNAs was assayed with cancer miRNA array respectively and detected with chemiluminescence imaging system. The internal control RNU48 was indicated in the image.

## Trouble Shooting

### Signals are too weak

- Total RNA may not contain small RNA
- RNA may be degraded
- If signals of the alignment spots are weak as well, the incubation of the membrane with Streptavidin HRP conjugate may be too short or the exposure time may be too short.

### Uneven background

- Substrate may not evenly overlaid on the membrane

## References

1. Ambros, V. (2003), MicroRNA pathways in flies and worms: growth, death, fat, stress, and timing. *Cell*. 113:673–676.
2. Poy, et. al. (2004) A pancreatic islet-specific microRNA regulates insulin secretion. *Nature*. 432:226–230.
3. Berezikov, et. al. (2006) Many novel mammalian microRNA candidates identified by extensive cloning and RAKE analysis. *Genome Res*. 16:1289–1298.
4. Calin, et. al. (2002) Frequent deletions and down-regulation of micro-RNA genes *miR15* and *miR16* at 13q14 in chronic lymphocytic leukemia. *Proc. Natl. Acad. Sci. USA*. 99:15524–15529.
5. Calin, et. al. (2004) MicroRNA profiling reveals distinct signatures in B cell chronic lymphocytic leukemias. *Proc. Natl. Acad. Sci. USA*. 101:11755–11760.

### Schematic diagram of cancer miRNA array

Let-7a	Let-7b	Let-7c	Let-7d	Let-7e	Let-7f	Let-7g	Let-7i	miR-1	miR-7
miR-9	miR-10a	miR-15a	miR-15b	miR-16	miR-17-5p	miR-18a	miR-18b	miR-19a	miR-19b
miR-20a	miR-21	miR-25	miR-28	miR-34a	miR-99a	miR-122a	miR-124a	miR-125a	miR-125b
miR-126	miR-131	miR-133a	miR-133b	miR-143	miR-145	miR-146a	miR-146b	miR-148a	miR-155
miR-181a	miR-181b	miR-181c	miR-182	miR-192	miR-194	miR-195	miR-199a	miR-199b	miR-199a*
miR-200a	miR-200c	miR-204	miR-206	miR-216	miR-223	miR-224	miR-342	miR-368	miR-375
miR-9-1	miR-10b	miR-17-3p	miR-22	miR-23a	miR-24	miR-26a	miR-26b	miR-27a	miR-27b
miR-29a	miR-29b	miR-29c	miR-30a	miR-30a	miR-30b	miR-30c	miR-92	miR-92b	miR-93
miR-95	miR-101-1	miR-103	miR-106a	miR-106b	miR-107	miR-128a	miR-128b	miR-132	miR-134
miR-135b	miR-136	miR-137	miR-140	miR-141	miR-142-3p	miR-149	miR-150	miR-151	miR-153
miR-154	miR-181d	miR-183	miR-185	miR-186	miR-188	miR-190	miR-191	miR-196a	miR-196b
miR-197	miR-198	miR-200b	miR-202	miR-203	miR-205	miR-210	miR-214	miR-215	miR-218
miR-219	miR-221	miR-222	miR-296	miR-372	miR-373	miR-488	miR-100	miR-127	miR-142-5p
miR-31	miR-213	RNU48							