



## Transcription Factor Western Blot Detection Assay Kit

Catalog Number WB-0001

(For Research Use Only)

### Introduction

Western Blotting has become one of the most widespread techniques for detecting a specific protein in a sample. Although this technique has become standard procedure for many labs, each new antibody requires a great deal of optimization. Signosis Transcription Factor Western Blot assay kits provide a pre-optimized, simple, and highly sensitive tool to conduct the assay. These kits contain all the components needed to probe and detect 5 whole blots.

### Materials provided with the kit

- 100 µl HRP-labeled Anti-Rabbit Secondary Antibody (4° C)
- 100 µl HRP-labeled Anti-Mouse Secondary Antibody (4° C)
- 40 ml of 2x Blocking buffer (4° C)
- 40 ml of 5x Diluting buffer (4° C)
- 40 ml 10x Washing buffer (RT)
- 5 Detection sheets
- 2.5 ml Substrate A (4° C)
- 2.5 ml Substrate B (4° C)

### Materials and equipment needed

1. Primary Antibody (WA-0XXX or WA-1XXX, Signosis)
2. Shaker
3. Imaging system or X-ray film

### Reagent preparation before experiment

- Dilute the 2x Blocking buffer, the 5x Diluting buffer, and the 10x Washing buffer to 1x buffer with ddH<sub>2</sub>O

### Assay Procedure

#### 1. Blotting

- (1) After performing sample preparation, gel electrophoresis, and transfer, block membrane in blocking buffer overnight at 4° C with shaking.
- (2) Rinse with washing buffer for 5 minutes with shaking.
- (3) Incubate with primary antibody (1:1000) in diluting buffer with shaking.
- (4) Rinse with washing buffer for 10 minutes, 5 minutes, and 5 minutes with shaking.
- (5) Incubate with secondary antibody (1:5000) in diluting

- buffer for 1 hour or more with shaking.
- (6) Rinse 3 times with washing buffer for 10 minutes each with shaking.
- (7) The membrane is ready for developing.

#### 2. Detection

- (1) Slightly dry the membrane by touching a corner of the membrane on a paper towel.
- (2) Mix 0.5 ml of Substrate A and 0.5 ml of Substrate B.
- (3) Place the membrane on the bottom side of detection sheet on a flat surface, and overlay the membrane with 1ml of substrate solution. Ensure that the substrate is evenly distributed over the membrane.
- (4) Gently place the top side of detection sheet over the membrane without trapping air bubbles in the membrane.
- (5) Incubate at room temperature for 5 minutes.
- (6) Remove excess substrate by gently applying pressure over the top sheet using a paper towel.
- (7) Expose the membranes using either Hyperfilm (2-10 min) or a chemiluminescence imaging system (i.e., FluorChem imager from Alpha Innotech). With either method, experiment with different exposure times.

### Sample of TF Western Blot analysis

