



## Mouse Obesity ELISA Strip for Profiling 8 Cytokines

Catalog Number EA-1071

(For Research Use Only)

### Introduction


Obesity increases the risk for the metabolic syndrome, diabetes, hypertension, atherosclerosis, and thrombosis (1-5). A number of proteins have been identified to be relevant to the development of the metabolic syndrome, diabetes, and cardiovascular disease with obesity. These include adiponectin, leptin, TNF $\alpha$ , IGF-1, Resistin, TGF $\beta$ , IL-6, and PAI-1. Plasma concentrations of these proteins are usually measured by ELISA. To systematically examine the effects, Signosis developed an ELISA strip which allows simultaneous determination of 8 proteins: leptin, TNF $\alpha$ , IGF-1, IL-6, VEGF, IL-1a, IL-1b, and MCP-1. The difference of these proteins between two samples can be determined through data comparison. Therefore, it facilitates the discovery of the change of these proteins in different samples.

### Principle of the assay

In each well of the strip, a primary antibody against a specific obesity cytokine is coated and 8 wells of the strip are coated with 8 different antibodies. Therefore, total 8 wells of a strip allow measurement of 8 different cytokines. The test sample is allowed to react simultaneously with pairs of two antibodies, resulting in the obesity cytokines being sandwiched between the solid phase and enzyme-linked antibodies. After incubation, the wells are washed to remove unbound-labeled antibodies. A HRP substrate, TMB, is added to result in the development of a blue color. The color development is then stopped with the addition of Stop Solution changing the color to yellow. The concentrations of obesity cytokines are directly proportional to the color intensity of the test sample. Absorbance is measured spectrophotometrically at 450 nm.

### Materials provided with the kit

- 12 strips, each coated with 8 different antibodies against mouse obesity cytokines (4°C).
- Biotin labeled antibody mixture against 8 different obesity cytokines (-20°C).
- Streptavidin-HRP conjugate (4°C)
- 1X Diluent buffer (4°C)
- 5X Assay wash buffer (RT)
- Substrate (4°C)
- Stop Solution (4°C)



Incubate with Detection antibody mixture

Incubate with HRP-Streptavidin

Add substrate TMB

OD450 reading




Diagram of Mouse Obesity ELISA Strip Analysis

### Reagent preparation before starting experiment

- Dilute the 5x Assay wash buffer to 1x buffer  
40ml 5x Assay wash buffer  
160ml ddH<sub>2</sub>O
- Use original or 10-fold diluted serum-free conditioned media, cell lysates, or sera. Samples can be diluted with 1 X Diluent buffer.
- Dilute 50 times of biotin labeled antibody mixture with 1X Diluent buffer.
- Dilute 200 times of streptavidin-HRP with 1X Diluent buffer.

## Recommendation

- The product intends to be used for comparison of 12 different samples. The differences of the obesity markers among the samples can be easily identified and determined.
- If you would like to quantitatively measure the proteins in the samples, please order EA-1072. It is protein standards which can be used for making standard curves through a series of 2-fold dilutions.


## Assay procedure

1. Cut the sealing film over the plate and remove it from the desired number of well strips. Make sure the rest of wells are well sealed.
2. Add 100µl of Standard, control, or sample per well and incubate for 1 hour at room temperature with gentle shaking.
3. Aspirate each well and wash by adding 200µl of 1X Assay wash buffer. Repeat the process three times for a total of three washes. Complete removal of liquid at each wash. After the last wash, remove any remaining liquid by inverting the plate against clean paper towels.
4. Add 100µl of diluted biotin-labeled antibody mixture to each well and incubate for 1 hour at room temperature with gentle shaking.
5. Repeat the aspiration/wash as in step 3.
6. Add 100µl of diluted streptavidin-HRP conjugate to each well and incubate for 45 min at room temperature with gentle shaking.
7. Repeat the aspiration/wash as in step 3.
8. Add 100µl substrate to each well and incubate for 5-30 minutes.
9. Add 50µl of Stop solution to each well. The color in the wells should change from blue to yellow.
10. Determine the optical density of each well with a microplate reader at 450 nm within 30 minutes.

## References

1. Moller DE, Flier JS. Insulin resistance—mechanisms, syndromes, and implications. *N Engl J Med* 1991;325:938–48.
2. Hubert HB, Feinleib M, McNamara PM, Castelli WP. Obesity as an independent risk factor for cardiovascular disease: a 26-year follow-up of participants in the Framingham Heart Study. *Circulation* 1983;67:968–77.
3. Loskutoff DJ, Samad F. The adipocyte and hemostatic balance in obesity. *Studies of PAI-1. Arterioscler Thromb Vasc Biol* 1998;18: 1–6.
4. Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults. Executive summary of the third report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). *JAMA* 2001;285:2486–97.
5. Ford ES, Giles WH, Dietz WH. Prevalence of the metabolic syndrome among US adults: findings from the third National Health and Nutrition Examination Survey. *JAMA* 2002;287:356–9.

## Example of standard curve



	Black bar	Grey bar
<b>Leptin</b>	2ng/ml	0.5ng/ml
<b>TNFa</b>	10ng/ml	2.5ng/ml
<b>IGF1</b>	10ng/ml	2.5ng/ml
<b>IL6</b>	4ng/ml	1ng/ml
<b>VEGF</b>	2ng/ml	0.5ng/ml
<b>IL-1a</b>	2ng/ml	0.5ng/ml
<b>IL-1b</b>	4ng/ml	1ng/ml
<b>MCP-1</b>	2ng/ml	0.5ng/ml