

## Malachite Green Phosphate Assay

Cat. No. 8118

2500 tests

### Introduction

---

Phosphorylation and dephosphorylation, which are the addition and removal of phosphate ( $\text{PO}_4^{3-}$ ) group to and from protein molecules respectively, are important regulatory mechanisms involved in cell-cycle regulation and signal transduction. ScienCell<sup>TM</sup> Malachite Green Phosphate Assay provides a simple colorimetric method for the determination of inorganic, soluble phosphate concentrations based on the complexation of malachite green oxalate with phosphomolybdate under acidic conditions. Applications of the kit include measurement of phosphate release, quantification of protein- or lipid- phosphorylation and etc. For quantification of protein or lipid bound phosphate, the phosphorylated proteins or lipids need to be hydrolyzed and neutralized before phosphate measurement are performed.

### Kit Components

Cat. No.	# of vials	Name	Quantity	Storage
8118a	1	10 mM Phosphate Standard	0.25 ml	4°C
8118b	1	Malachite Green Reagent A	25 ml	4°C
8118c	1	Malachite Green Reagent B	25 ml	4°C

### Quality Control

---

Data from ScienCell<sup>TM</sup> Malachite Green Phosphate Assay of phosphate solutions with concentrations ranging from 1.5625 to 37.5  $\mu\text{M}$  shows a linear relationship between  $\text{OD}_{630\text{nm}}$  and phosphate concentration (Figure 1).

### Procedures

---

#### A. Preparation of phosphate standards:

1. Prepare a 50  $\mu\text{M}$  phosphate standard by adding 5  $\mu\text{l}$  of 10 mM Phosphate Standard to 995  $\mu\text{l}$  of DI  $\text{H}_2\text{O}$ .
2. Prepare a phosphate standard curve using the serial dilutions of the 50  $\mu\text{M}$  phosphate standard according to Table 1. 160  $\mu\text{l}$  of phosphate solution is prepared for each point to provide three replicates of 50  $\mu\text{l}$ .

#### B. Assay procedure:

1. Add 50  $\mu\text{l}$  of phosphate standard, sample or blank into each well of 96 well plate.
2. Add 10  $\mu\text{l}$  of Malachite Green Reagent A to each well, mix and incubate for 10 minutes at room temperature.
3. Add 10  $\mu\text{l}$  of Malachite Green Reagent B to each well, mix and incubate for 20 minutes at room temperature.
4. Read absorbance at 630 nm on a plate reader.

### C. Calculations:

1. Average the OD<sub>630nm</sub> of replicate wells of each phosphate standard, sample and blank. Subtract the average OD<sub>630nm</sub> value of the blank from the average OD<sub>630nm</sub> values obtained with all other samples.
2. Based on the calibrated OD<sub>630nm</sub> of the phosphate standard, make a standard curve by plotting OD<sub>630nm</sub> as a function of phosphate concentration. (See Figure 1 for a typical standard curve.) Determine the equation and R<sup>2</sup> value of the trend line.
3. Suppose the equation of the trend line of the standard curve is  $y = Ax + B$ , calculate the phosphate concentration of samples as follows:

$$[PO_4^{3-}] = \frac{OD_{630nm} - B}{A}$$

No.	50 μM phosphate (μL)	DI H <sub>2</sub> O (μL)	Phosphate concentration (μM)
1	160	0	50
2	80	80	25
3	40	120	12.5
4	20	140	6.25
5	10	150	3.125
6	5	155	1.5625
7	0	160	0 (Blank)

Table 1. Preparation of phosphate standards in Malachite Green Phosphate Assay.




Figure 1. A typical phosphate standard curve measured by ScienCell™ Malachite Green Phosphate Assay.