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(K4301, K4302, K4303, K4304)

Description

Lumi-Glow™ luminescent kinase assay system is the homogeneous ready to use reagent for quantitatively detecting kinase levels in the samples. It featured high signal to noise ratio, good reproducibility and stability. The signal half-time is over 2 hours. The homogeneous formulation eliminates the procedures of addition of the reagent after cell lysis, as well as decreases experiment errors generated by multiple additions. Moreover, the product is suitable for high throughput sample measurement due to the stable “glow” type luminescence.

Reagent components

Luciferase, luciferin, buffer were mixed and filled in the 10 ml or 100 ml reagent bottles, three packages are showed below:

Cat#	Size	Assays of 96-well plate	Assays of 384-well plate
K4301	10 ml	200 wells	1,000 wells
K4302	100 ml	2,000 wells	10,000 wells
K4303	10X100 ml	20,000 wells	100,000 wells

Procedures

Optimization of kinase reaction

To keep Lumi-Glow™ luminescent kinase assay system in better conditions, we recommend kinase reaction conditions should be optimized, including concentrations of ATP, kinases, and substrates. For 96-well plate, we recommend add 50 µl of Lumi-Glow™ luminescent kinase reagent into 50 µl of kinase reaction system; for 384-well plate, 10 µl of Lumi-Glow™ luminescent kinase reagent into 10 µl of kinase reaction system; for other plates, add 1 volume of Lumi-Glow™ luminescent kinase reagent into 1 volume of kinase reaction system

Optimization of ATP concentration

- Dilute ATP with kinase reaction buffer at 1:2 ratio, add into the plate, add sufficient kinase and kinase substrate. For control wells, add same concentrations of ATP without kinase or kinase substrate. Mix and incubate at room temperature for appropriate time
- Add equal volume of Lumi-Glow™ luminescent kinase assay system, mix
- Read the plate in “luminescence” plate reader after 10 min of reaction

Optimization of kinase substrate concentration

- Dilute kinase substrate with kinase reaction buffer at 1:2 ratio, add into the plate, add sufficient kinase and optimal ATP. Mix and incubate at room temperature for appropriate time
- Add equal volume of Lumi-Glow™ luminescent kinase assay system, mix
- Read the plate in “luminescence” plate reader after 10 min of reaction

Optimization of kinase concentration

- Dilute kinase with kinase reaction buffer at 1:2 ratio, add into the plate, add optimal ATP and kinase substrate. Mix and incubate at room temperature for appropriate time
- Add equal volume of Lumi-Glow™ luminescent kinase assay system, mix
- Read the plate in “luminescence” plate reader after 10 min of reaction

Protocols (for 96-well plate)

- Add ATP, kinase, kinase substrate, kinase reaction buffer into the 96-well plate. Total volume is 50 µl, mix, incubate at room temperature for appropriate time
- Add 50 µl of Lumi-Glow™ luminescent kinase assay system, mix
- Read the plate in “luminescence” plate reader after 10 min of reaction

Storage conditions

Store Lumi-Glow™ luminescent kinase assay kit at -20°C. The product can be stored at 4°C for one week, and can tolerate 10 times of “freeze and thaw”