Lumi-Glow™ *Renilla* Luciferase Assay System

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(F2101, F2102, F2103)

Description

Lumi-Glow™ *Renilla* luciferase assay system is the easy to use reagent for quantitatively detecting *Renilla* luciferase levels in the cells or other samples. It featured high signal to noise ratio, good reproducibility and stability. The signal half-time is over 2 hours. The convenient formulation eliminates the procedures of addition of the reagent after cell lysis, as well as decreases experiment errors generated by multiple additions. Moreover, the product is suitable for high throughput sample measurement due to the stable “glow” type luminescence.

Reagent components

Buffer was filled in the 10 ml or 100 ml reagent bottles, and substrate was filled in 2 ml tubes, three packages are showed below:

Cat#	Size	Assays for 96-well plate	Assays for 384-well plate
R6501	Buffer, 10 ml	200 wells	800 wells
	Substrate, 100 µl		
R6502	Buffer, 100 ml	2,000 wells	8,000 wells
	Substrate, 1 ml		
R6503	Buffer, 10X100 ml	20,000 wells	80,000 wells
	Substrate, 10x1 ml		

Procedures

- Reagent preparation: add 1 volume of Lumi-Glow™ *Renilla* substrate to 100 volume of Lumi-Glow™ *Renilla* buffer, i.e. add 50 µl of substrate to 5 ml of buffer to reconstitute the detection reagent
- Equilibrate the reagent and samples at RT for 10-60 min till the reagent temperature reaches room temperature. Gently mix the reagent through inverting the bottle for several times
- Add 50 µl of reagent to 100 µl of cells or samples in 96-well plate or add 10 µl of reagent to 20 µl of cells or samples in 384-well plate
- Gently mix the cells or samples
- Read the plate in the appropriate luminescence plate reader after 10 min of reagent addition

Storage conditions

Store Lumi-Glow™ *Renilla* luciferase assay kit at -20°C. It is strongly recommended that substrate and buffer are freshly mixed according to the experiment needs