



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

INSTRUCTION MANUAL

Quick-RNATM MicroPrep

Catalog Nos. **R1050 & R1051**

Highlights

- High-quality total RNA from a wide range of samples – single to 10⁶ cells.
- Isolate small and large RNAs into separate fractions (optional).
- *DNA-free* RNA for use in any downstream application. *DNase I included.*

Contents

Product Contents	1
Specifications.....	1
Product Description.....	2
Buffer Preparation	3
Protocols.....	3, 4
Appendices	5
Ordering Information	6
Related Products	7

Satisfaction of all Zymo Research products is guaranteed. If you should be dissatisfied with this product please contact us.

Some difficult-to-lyse samples may require mechanical or enzymatic homogenization. For assistance, contact us at tech@zymoresearch.com.

For 10^2 to 10^7 cells, use the **Quick-RNA™ MiniPrep** (Cat. Nos. R1054, R1055).

Product Contents

Quick-RNA™ MicroPrep (Kit Size)	R1050 (50 Preps.)	R1051 (200 Preps.)
RNA Lysis Buffer	50 ml	2x 100 ml
RNA Prep Buffer	25 ml	100 ml
RNA Wash Buffer¹ (concentrate)	24 ml	2x 48 ml
DNase/RNase-Free Water	6 ml	2x 10 ml
DNase I Set² DNase I (250 U) & 10x DNase I Reaction Buffer (1 ml)	1 set	4 sets
Zymo-Spin™ IC Columns	50	200
Collection Tubes	50	200
Instruction Manual	1	1

Note - Integrity of kit components is guaranteed for up to one year from date of purchase. Reagents are routinely tested on a lot-to-lot basis to ensure they provide the highest performance and reliability.

Storage Temperature - Store all kit components (*i.e.*, buffers, columns) at room temperature. Store reconstituted DNase I at -20 °C.

¹ Add 96 ml 100% ethanol (104 ml 95% ethanol) to the 24 ml **RNA Wash Buffer** concentrate or 192 ml 100% ethanol (208 ml 95% ethanol) to the 48 ml **RNA Wash Buffer** concentrate before use.

² Add 275 µl **DNase/RNase-Free Water** per vial to reconstitute the lyophilized **DNase I** (E1009) at 1 U/µl. Mix by gentle inversion. Store aliquots at -20 °C.

Specifications

- **Sample Sources** – Cells or tissue samples, yeast, plant, bacteria, buccal cells, buffy coat, plasma, serum, and other biological liquids. *Compatible with RNA Shield™ (Cat. No. R1100) and RNAlater™.*
- **Sample Storage** – Samples homogenized in RNA Lysis Buffer are stable and can be stored frozen prior to purification.
- **Sample Size** – Up to 10^6 cells or 5 mg tissue.
- **RNA Purity** – High quality RNA ($A_{260}/A_{280} > 1.8$, $A_{260}/A_{230} > 1.8$) suitable for all downstream RNA-based manipulations.
- **RNA Recovery** – Up to 10 µg RNA can be eluted into ≥6 µl RNase-free water allowing for a highly concentrated sample.
- **RNA Storage** – RNA is eluted with RNase-free water and can be stored frozen. RNase inhibitors can be included for prolonged storage.
- **Equipment Needed** – Microcentrifuge.

Note - [™] Trademarks of Zymo Research Corporation. This product is for research use only and should only be used by trained professionals. It is not intended for use in diagnostic procedures. Some reagents included with this kit are irritants. Wear protective gloves and eye protection. Follow the safety guidelines and rules enacted by your research institution or facility. *RNAlater™* is a trademark of Ambion, Inc., Austin, Texas and is protected by various U.S. and foreign patents.


ZYMO RESEARCH CORP.

Phone: (949) 679-1190 • Toll Free: (888) 882-9682 • Fax: (949) 266-9452 • info@zymoresearch.com • www.zymoresearch.com


Product Description

The **Quick-RNA™ MicroPrep** kit is an innovative product designed for the easy, reliable, and rapid isolation of DNA-free RNA from a wide range of cell (*up to 10⁶*) and tissue samples (*up to 5 mg*). The procedure combines a unique buffer system with *Fast-Spin* column technology to yield high quality total RNA (*including small RNAs 17-200 nt*) in about 10 minutes.

The procedure is simple: Add the provided **RNA Lysis Buffer** to a sample, then purify the RNA using the **Zymo-Spin™ Columns**. The result is highly-concentrated, *DNA-free* RNA that is suitable for subsequent RNA-based methods including RT-PCR, hybridization, sequencing *etc.* In addition, the kit can be used for enrichment of small and large RNAs in two separate fractions (see Appendix, page 5).



The **Quick-RNA™ MicroPrep** yields high quality total RNA. High levels of genomic DNA contamination are present in the preps from Suppliers Q & P but not with the **Quick-RNA™ MicroPrep**. Total RNA was isolated from human epithelial cells (sans DNase treatment).




RNA isolated with the **Quick-RNA™ MicroPrep** is DNA-free. Samples isolated with Supplier Q's kit are provided for comparison. Total RNA was isolated from 10⁶ human epithelial cells (with in-column DNase treatments for both kits). Each amplification curve represents an average of three independent isolation experiments.

For **Assistance**, please contact Zymo Research Technical Support at 1-888-882-9682 or e-mail tech@zymoresearch.com.

Notes:

Use the **Direct-zol™ RNA MiniPrep** (Cat. Nos. R2050, R2051, R2052, R2053) for isolation of RNA *directly* (without phase separation) from samples in Trizol®, *etc.*

Use the **RNA Shield™** (Cat. Nos. R1100-50, R1100-250) for safe sample storage and transport at ambient temperatures.



The **Quick-RNA™** kits yield high quality RNA with high "RNA Integrity Numbers" (2200 TapeStation, Agilent).

ZYMO RESEARCH CORP.

Phone: (949) 679-1190 • Toll Free: (888) 882-9682 • Fax: (949) 266-9452 • info@zymoresearch.com • www.zymoresearch.com

Ensure the RNA isolation procedure is performed in an RNase-free environment.

Notes:

Samples homogenized in **RNA Lysis Buffer** can be stored frozen for processing at a later time.

ZR Bashing Bead™ Lysis Tubes are available separately (Cat. Nos. S6002, S6003).

Processing plant tissue and other samples containing polyphenolics, humic acids, melanin, etc. may require use of the **OneStep™ PCR Inhibitor Removal Kit** (Cat. No. D6030).

For whole-blood samples, use the **ZR Whole-Blood RNA MiniPrep™** (Cat. Nos. R1020, R1021).

Buffer Preparation

- ✓ Before starting, add 96 ml 100% ethanol (104 ml 95% ethanol) to the 24 ml **RNA Wash Buffer** concentrate (R1050) or 192 ml 100% ethanol (208 ml 95% ethanol) to the 48 ml **RNA Wash Buffer** concentrate (R1051).
- ✓ Add 275 µl **DNase/RNase-Free Water** per vial to reconstitute the lyophilized **DNase I** (E1009) at 1 U/µl. Mix by gentle inversion. Store aliquots at -20 °C.

Protocols

The RNA isolation consists of three steps: (I) *Sample Lysis/Homogenization*, (II) *Sample Clearing* and (III) *RNA Purification*.

All steps should be performed at room temperature (20-30 °C).

I. Sample Lysis/Homogenization

Recommended **RNA Lysis Buffer** volumes

RNA Lysis Buffer	100 µl	300 µl
Cells	Up to 10 ⁵	Up to 10 ⁶
Tissue	-	Up to 5 mg

Adherent Cells

Lyse cells directly in the culture container by removing liquid medium and adding **RNA Lysis Buffer** directly to the monolayer. Remove cells from culture surface by pipetting, scraping, etc.

Cells in Suspension

Pellet cells (≤500 x g), remove the supernatant completely then resuspend the cell pellet in **RNA Lysis Buffer**. Vortex briefly.

Tissue and Tough-to-Lyse Samples

Fresh or frozen tissue (animal, plant, insect, yeast or bacteria) can be mechanically homogenized (**ZR BashingBead™ Lysis Tubes**, Dounce or similar) directly in the **RNA Lysis Buffer**. Alternatively, tough-to-lyse samples can undergo enzymatic treatment (Proteinase K (Appendix C, page 5) – tissue, Zymolyase - yeast) prior to adding the **RNA Lysis Buffer**.

Liquids

DNase-treated RNA, labeling reactions, aqueous phase (following Trizol® extraction), and biological liquids can be processed directly by adding 4 volumes of **RNA Lysis Buffer** to each volume of sample (4:1) then mixing well. Plasma and serum samples require digestion with Proteinase K (Appendix C, page 5) before processing.

Samples in RNA Shield™

Bring samples homogenized and stored in **RNA Shield™** to room temperature (20-30 °C) before proceeding with *Sample Clearing* step.

Samples in RNA/ater™

To process cells or liquids in **RNA/ater™** (without reagent removal): Add 1 volume of RNase-free water or PBS to the sample (1:1). Then add 4 volumes **RNA Lysis Buffer** (4:1) and mix.

Note: Alternatively, remove the **RNA/ater™** then proceed with the recommended *Sample Lysis/Homogenization* (above) according to the sample type.

ZYMO RESEARCH CORP.

Phone: (949) 679-1190 • Toll Free: (888) 882-9682 • Fax: (949) 266-9452 • info@zymoresearch.com • www.zymoresearch.com

II. Sample Clearing

The following is recommended for cells and tissue (animal/plant) but can be omitted for cell-free liquids and low input samples ($\leq 10^5$ cells).

For particulate removal, centrifuge lysates at $\geq 12,000 \times g$ for 1 minute. Then transfer the supernatant into an RNase-free tube (*not provided*).

III. RNA Purification

All centrifugation steps should be performed at 10,000-16,000 $\times g$.

1. Add 1 volume ethanol (95-100%) to the sample in **RNA Lysis Buffer** (1:1). Mix well.
2. Transfer the mixture to a **Zymo-Spin™ IC Column**¹ in a **Collection Tube** and centrifuge for 30 seconds. Discard the flow-through.
3. **In-column DNase I Treatment** (optional)

This step can be used for trace DNA removal.

- a. Add 400 μ l **RNA Wash Buffer** to the column and centrifuge for 30 seconds. Discard the flow-through.
- b. Add 50 μ l of the **DNase I reaction mix** (see below) directly to the column matrix. Incubate the column at room temperature (20-30 °C) for 15 minutes, then centrifuge for 30 seconds.

For each sample to be treated, prepare **DNase I reaction mix**² in an RNase-free tube (not provided). Mix gently.

DNase I	5 μ l (1 U/ μ l)
10X DNase I Reaction Buffer	5 μ l
RNA Wash Buffer (with ethanol added)	40 μ l

4. Add 400 μ l **RNA Prep Buffer** to the column and centrifuge for 30 seconds. Discard the flow-through.
5. Add 700 μ l **RNA Wash Buffer** to the column and centrifuge for 30 seconds. Discard the flow-through.
6. Add 400 μ l **RNA Wash Buffer** and centrifuge the column for 2 minutes to ensure complete removal of the wash buffer.
7. Place the column into an RNase-free tube (*not provided*). Add ≥ 6 μ l **DNase/RNase-Free Water**³ directly to the column matrix, then centrifuge at top speed for 30 seconds.

Eluted RNA can be used immediately or stored frozen.

Notes:

¹ To process samples >800 μ l, **Zymo-Spin™** columns may be reloaded.

² When adjusting the volume and composition, make sure the **RNA Wash Buffer** in the **DNase I reaction mix** remains at 80% (v/v).

Unit definition - one unit increases the absorbance of a high molecular weight DNA solution at a rate of 0.001 A_{260} units/min/ml of reaction mixture at 25°C.

³ To maximize RNA yield, preheat the DNase/RNase-Free Water to 95° C, increase the elution volume and/or repeat the elution.

ZYMO RESEARCH CORP.

Phone: (949) 679-1190 • Toll Free: (888) 882-9682 • Fax: (949) 266-9452 • info@zymoresearch.com • www.zymoresearch.com

Notes:

¹ Minimal recommended sample volume is 50 µl. Adjust samples <50 µl with RNase-free water.

² To process samples >800 µl, **Zymo-Spin™** columns may be reloaded.

³ One unit (1 U) of enzyme will hydrolyze urea-denatured hemoglobin to produce 1.0 µmole of tyrosine per minute at pH 7.5 at 37°C.

Appendices

A. Purification of Small and Large RNAs into Separate Fractions

The **Quick-RNA™ MicroPrep** also allows for purification of small (17-200 nt) and large RNAs (>200 nt) into separate fractions. This procedure is compatible only with cell inputs (up to 10⁶) or previously isolated RNA. Processing requires two Zymo-Spin™ IC columns per prep.


All centrifugation steps should be performed at 10,000-16,000 x g.

1. Prepare a modified **RNA Lysis Buffer** with the addition of 1 volume ethanol (95-100%) (1:1). Mix. Add 2 volumes of this modified buffer to each volume of RNA sample¹ (2:1) or 300 µl buffer to pelleted animal cells. Mix well. Do not centrifuge the sample.
2. Transfer the mixture to the **Zymo-Spin™ IC Column** in a **Collection Tube** and centrifuge for 30 seconds. ***Important: Save the flow-through!***

To purify large RNAs (>200 nt) from this column, continue to *RNA Purification* (page 4, step 3).

3. To the **flow-through**, add 1 volume ethanol (95-100%) (1:1) and mix well.
4. Transfer the mixture to a **new Zymo-Spin™ IC Column** in a **Collection Tube** and centrifuge for 30 seconds. Discard the flow-through.

To purify small RNAs (17-200 nt) from this column, continue to *RNA Purification* (page 4, step 3)².



Total RNA (>17 nt), large (>200 nt) or small RNAs (17-200 nt) are effectively partitioned and purified with the **Quick-RNA™** kit.

B. RNA Clean-Up

The **Quick-RNA™ MicroPrep** can also be used for RNA clean-up. To clean RNA following DNase treatment, labeling reactions or from an aqueous phase (following acid-guanidinium-phenol extractions); simply add 4 volumes **RNA Lysis Buffer** and proceed with the *RNA Purification* (page 4).

C. Proteinase K Digestion

Tissue, Plasma, Serum
(RNase-free water for solid tissue up to 5 mg)
2X Digestion Buffer
Proteinase K

95 µl

95 µl

≥6 U³

Prepare a Proteinase K reaction mix (see example above, scale-up as necessary). Incubate at 55°C for 30 minutes (plasma, serum) or 1-3 hours (tissue). Then add 4 volumes **RNA Lysis Buffer** and proceed to *Sample Clearing* (page 4).

ZYMO RESEARCH CORP.

Phone: (949) 679-1190 ▪ Toll Free: (888) 882-9682 ▪ Fax: (949) 266-9452 ▪ info@zymoresearch.com ▪ www.zymoresearch.com

Ordering Information

Product Description	Input	Binding	Catalog No.	Kit Size
Quick-RNA™ MicroPrep	~1-10 ⁶ cells	~10 µg	R1050 R1051	50 Preps. 200 Preps.
Quick-RNA™ MiniPrep	~10 ² -10 ⁷ cells	~100 µg	R1054 R1055	50 Preps. 200 Preps.
Quick-RNA™ MidiPrep	~10 ⁶ -10 ⁸ cells	~1 mg	R1056	25 Preps.
ZR-96 Quick-RNA™	~1-10 ⁶ cells	~10 µg/well	R1052 R1053	2x 96 Preps. 4x 96 Preps.

For Individual Sale	Catalog No.	Amount
RNA Lysis Buffer	R1060-1-50 R1060-1-100	50 ml 100 ml
RNA Prep Buffer	R1060-2-10 R1060-2-25 R1060-2-100	10 ml 25 ml 100 ml
RNA Wash Buffer (concentrate)	R1003-3-6 R1003-3-12 R1003-3-24 R1003-3-48	6 ml 12 ml 24 ml 48 ml
DNase I (lyophilized) (250 U supplied with 10x DNase I Reaction Buffer)	E1009	1 set
Zymo-Spin™ IC Column	C1004-50-G C1004-250-G	50 250
Collection Tube	C1001-50 C1001-500 C1001-1000	50 500 1000
DNase/RNase-Free Water	W1001-1 W1001-6 W1001-10	1 ml 6 ml 10 ml

ZYMO RESEARCH CORP.

Phone: (949) 679-1190 ▪ Toll Free: (888) 882-9682 ▪ Fax: (949) 266-9452 ▪ info@zymoresearch.com ▪ www.zymoresearch.com

Related Products

Product	Description	Prep/Format	Catalog
RNA Clean-Up			
RNA Clean & Concentrator™ -5	Cleanup and concentration of modified, labeled, impure, diluted, DNase treated RNA (≥ 17nt) and purification of RNA from aqueous phase of organic extracts. <i>Note: DNA-free RNA Kit™ includes DNase I</i>	50/column 200/column	R1015 R1016
RNA Clean & Concentrator™ -25		50/column 100/column	R1017 R1018
RNA Clean & Concentrator™ -100		25/column	R1019
ZR-96 RNA Clean & Concentrator™		2x 96/plate	R1080
DNA-Free RNA Kit™		50/column 200/column	R1013 R1014
Oligo Clean & Concentrator™	Cleanup and concentration of RNA and/or DNA oligos. Good for clean-up of miRNAs and siRNAs.	50/column 200/column	D4060 D4061
ZR-96 Oligo Clean & Concentrator™		2x 96/plate 4x 96/plate	D4062 D4063
ssDNA/RNA Clean & Concentrator™	Separation of short ssRNA and ssDNA (up to 200 nt) from double stranded species.	20/column 50/column	D7010 D7011
Zymoclean™ Gel RNA Recovery Kit	Recovery of RNA from agarose gels.	50/column	R1011
ZR small-RNA™ PAGE Recovery Kit	Small RNA (> 17nt) from polyacrylamide gels.	20/column	R1070
OneStep™ PCR Inhibitor Removal Kit	Removal of polyphenolics, humic/fulvic acids, tannins, melanin etc. from RNA.	50/column	D6030
OneStep™ -96 PCR Inhibitor Removal Kit		2x 96/plate	D6035
RNA from Samples in TRI Reagent (Small RNA Recovery)			
Direct-zol™ RNA MiniPrep	RNA (>17 nt) from TRI Reagent®, TRIzol®, and all other acid-guanidinium-phenol based reagents without phase separation. DNase I included.	50/column 200/column	R2050 R2052
Direct-zol™ RNA MiniPrep w/ TRI Reagent		50/column 200/column	R2051 R2053
Direct-zol™ -96 RNA		2x 96/plate 4x 96/plate	R2054 R2056
Direct-zol™ -96 RNA w/ TRI Reagent		2x 96/plate 4x 96/plate	R2055 R2057
Direct-zol™ -96 MagBead RNA	RNA (>17 nt) from TRI Reagent®, TRIzol®, and all other acid-guanidinium-phenol based reagents without phase separation. These kits are in a magnetic bead format that is adaptable for high-throughput and automated protocols. DNase I included.	2x 96/plate 4x 96/plate 8x 96/plate	R2100 R2102 R2104
Direct-zol™ -96 MagBead RNA w/ TRI Reagent		2x 96/plate 4x 96/plate 8x 96/plate	R2101 R2103 R2105
RNA from Cells & Tissue			
Quick-RNA™ MicroPrep	Total RNA (>17 nt) from cultured cells, tissue samples, buccal cells, buffy coat, plasma, serum, and other biological liquids. DNase I included.	50/column 200/column	R1050 R1051
Quick-RNA™ MiniPrep		50/column 200/column	R1054 R1055
Quick-RNA™ MidiPrep		25/column	R1056
ZR-96 Quick-RNA™		2x 96/plate 4x 96/plate	R1052 R1053
ZR-Duet™ DNA/RNA MiniPrep	Parallel purification of DNA/RNA from cells.	50/column	D7001
Pinpoint™ Slide RNA Isolation System Kit I	RNA from fresh/frozen tissue sections.	50/column	R1003
Pinpoint™ Slide RNA Isolation System Kit II	RNA from paraffin-embedded (FFPE) tissue.	50/column	R1007
RNA from Biological Liquids			
ZR Viral RNA Kit™	RNA (DNA) from body fluids (plasma, serum, CSF, urine).	50/column 200/column	R1034 R1035
ZR-96 Viral RNA Kit™		2x 96/plate 4x 96/plate	R1040 R1041
ZR Viral DNA/RNA Kit™		25/column 100/column	D7020 D7021
ZR Whole-Blood RNA MiniPrep™	RNA from whole blood or partitioned blood.	50/column 100/column	R1020 R1021
ZR Urine RNA Isolation Kit™	Cellular and endosomal RNA from urine.	20/column 50/column	R1038 R1039
RNA from Tough-to-Lyse Samples			
ZR Fungal/Bacterial RNA MicroPrep™	RNA from bacteria, yeast, fungi; BashingBead™ lysis.	50/column	R2010
ZR Fungal/Bacterial RNA MiniPrep™		50/column	R2014
ZR Plant RNA MiniPrep™	RNA from leaves, stems, buds, flowers, fruits, seeds, etc; BashingBead™ lysis, RT/PCR inhibitor removal.	50/column	R2024
ZR Tissue & Insect RNA MicroPrep™	RNA from insect, arthropod specimen and small tissue samples; BashingBead™ lysis.	50/column	R2030
ZR Soil/Fecal RNA MicroPrep™	RNA from soil, sludge, sediment, feces.	50/column	R2040
YeaStar RNA Kit™	RNA from yeast strains susceptible to Zymolyase.	50/column	R1002
RNA Sample Preservation and Storage			
RNA Shield™	Cells, biological liquid, tissue storage and RNA purification.	50 ml 250 ml	R1100-50 R1100-250
RNA Shield™ Purification Kit (RNA Shield™ reagent included)		50/column	R1100
RNA Shield™ Purification Kit (RNA Shield™ reagent is <u>not</u> included)		50/column	R1101
Enzymes and Markers			
DNase I w/ 10X Reaction Buffer	Lyophilized	250 U	E1009
ZR small-RNA™ Ladder	ssRNA (17, 21, 25, 29 nt)	10 µg	R1090

ZYMO RESEARCH CORP.

Phone: (949) 679-1190 ▪ Toll Free: (888) 882-9682 ▪ Fax: (949) 266-9452 ▪ info@zymoresearch.com ▪ www.zymoresearch.com