

# Yeast $\beta$ -Galactosidase Assay Kit

75768

0856.3

**Number****Description**

75768

**Yeast  $\beta$ -Galactosidase Assay Kit****Kit Contents:****Y-PER™ Yeast Protein Extraction Reagent**, 25mL, store at room temperature**2X  $\beta$ -Galactosidase Assay Buffer**, 25mL, store at -20°C **$\beta$ -Galactosidase Assay Stop Solution**, 25mL, contains 1 M Na<sub>2</sub>CO<sub>3</sub>, store at room temperature**Storage:** Upon receipt store individual components as indicated. Product is shipped with dry ice.**Introduction**

The Thermo Scientific Yeast  $\beta$ -Galactosidase Assay Kit allows for qualitative or quantitative determination of  $\beta$ -galactosidase activity in solution directly from colonies growing on solid medium. A portion of the colony is suspended in a mixture of Thermo Scientific Y-PER Yeast Protein Extraction Reagent and  $\beta$ -Galactosidase Assay Buffer. After a brief incubation period, the solution turns yellow from the hydrolysis of *o*-nitrophenyl- $\beta$ -D-galactopyranoside (ONPG) to *o*-nitrophenol (ONP) and galactose in a mildly alkaline solution. The assay becomes quantitative if cell quantity is first determined by measuring the optical density at 660nm.

The gene encoding  $\beta$ -galactosidase (*lacZ*) of *E. coli* has been widely used as a reporter gene in many different prokaryotic and eukaryotic organisms.  $\beta$ -Galactosidase activity is also commonly used as an indicator of protein:protein interactions *in vivo* using two-hybrid systems. The interaction strength is verified or quantitated using a  $\beta$ -galactosidase activity assay.

Common protocols for determining  $\beta$ -galactosidase activity from liquid cultures involve harvesting cells by centrifugation, washing cells several times, and lysing cells by either subjecting them to freeze/thaw cycles or by treating them with detergents or organic solvents. These traditional protocols are not amenable to projects that involve screening cells growing in 96-well plates. The Yeast  $\beta$ -Galactosidase Assay Kit allows cells obtained from solid media or liquid culture to be assayed directly with no harvesting or washing steps.

**Yeast Colony  $\beta$ -galactosidase Assay Protocol - Qualitative Method**

1. Thaw the 2X  $\beta$ -Galactosidase Assay Buffer on ice.
2. Add a volume of the assay buffer to an equal volume of Y-PER Reagent to make the working solution (WS). Each colony requires 50-100 $\mu$ L of the WS.
3. Pipette 50-100 $\mu$ L of the WS into each microcentrifuge tube. Prepare one tube for each colony to be assayed.
4. Use a sterile inoculation loop, toothpick or pipette tip to transfer a portion of a single colony to one of the microcentrifuge tubes containing the WS. Mix gently with a vortex mixer to create a homogeneous solution.
5. Incubate reaction tube at room temperature or 37°C (optimal) until a color change is observed.

**Note:** The solution will become yellow within minutes depending on the amount of  $\beta$ -galactosidase in the sample. If the color change is not apparent, measure the absorbance at 420nm against a blank containing the WS, which is useful for samples with low levels of  $\beta$ -galactosidase activity.

---

## Yeast Colony $\beta$ -galactosidase Assay Protocol - Quantitative Method

1. Thaw 2X  $\beta$ -Galactosidase Assay Buffer on ice.
2. Pipette 250 $\mu$ L of Y-PER Reagent into a microcentrifuge tube. Prepare one tube for each colony.
3. Use a sterile inoculation loop, toothpick or pipette tip to transfer a portion of the colony to one of the tubes containing the Y-PER Reagent. Mix gently with a vortex mixer to create a homogeneous solution.  
**Note:** For calculations, use 0.25mL for cell volume.
4. Determine the OD<sub>660</sub> of the solution and record this value.
5. Use a timer to monitor the reaction. Add 250 $\mu$ L of 2X  $\beta$ -Galactosidase Assay Buffer to the microcentrifuge tube and start timer.
6. Incubate the reaction tube at room temperature or 37°C (optimal) until a color change is observed.  
**Note:** The solution will become yellow within minutes depending on the amount of  $\beta$ -galactosidase in the sample.
7. When the color change appears, add 200 $\mu$ L of  $\beta$ -Galactosidase Assay Stop Solution to the reaction tube and vortex for 15 seconds. Stop the timer and record the total reaction time.
8. Remove the cell debris from the reaction tube by centrifuging at 13,000  $\times$  g for 30 seconds.
9. Transfer supernatant to a cuvette and measure the absorbance at 420nm against a blank containing 250 $\mu$ L of Y-PER Reagent, 250 $\mu$ L of the 2X  $\beta$ -Galactosidase Assay Buffer, and 200 $\mu$ L of the  $\beta$ -Galactosidase Assay Stop Solution.  
**Note:** The reaction time will vary depending on the level of  $\beta$ -galactosidase expression. Absorbance values between 0.02-1.0 are within the linear range of the assay.
10. Proceed to Calculations section.

## Yeast $\beta$ -galactosidase Microplate Plate Assay Protocol (non-stopped)

1. Thaw 2X  $\beta$ -Galactosidase Assay Buffer on ice.
2. Add a volume of this assay buffer with an equal volume of Y-PER Reagent to make the working solution (WS). Each well requires 100 $\mu$ L of the WS. The maximum culture volume allowed per well is 100 $\mu$ L.
3. Determine the OD<sub>660</sub> of all test cultures and record values. Include one well that contains growth media only (no cells) as a blank for the spectrophotometer.
4. Add 100 $\mu$ L of each culture to the wells.
5. Use a timer to monitor the reaction. Apply 100 $\mu$ L of the WS to each well and start timer.
6. Incubate the plate at room temperature for approximately 30 minutes.
7. Use the well containing medium only (no cells) to zero the instrument. Measure the absorbance at 420nm of each well.  
**Note:** The reaction time will vary depending on the level of  $\beta$ -galactosidase expression. Absorbance values between 0.02-1.0 are within the linear range of the assay.
8. Proceed to Calculations Section.

## Yeast $\beta$ -galactosidase Microplate Plate Assay Protocol (stopped)

1. Thaw 2X  $\beta$ -Galactosidase Assay Buffer on ice.
2. Add a volume of the assay buffer to an equal volume of Y-PER Reagent to make the working solution (WS). Each well requires 70 $\mu$ L of the WS. The maximum cell culture volume allowed per well is 70 $\mu$ L.
3. Determine the OD<sub>660</sub> of all test cultures and record values. Include one well that contains only the growth medium (no cells) as a blank for the instrument.
4. Add 70 $\mu$ L of each culture to individual wells in the microplate.
5. Use a timer to monitor the reaction. Apply 70 $\mu$ L of the WS to each well and start the timer.

6. When the yellow color appears, add 56µL of β-Galactosidase Assay Stop Solution to each well and mix for 15 seconds. Stop the timer and record the total reaction time.
7. Use the well containing medium only (no cells) to zero the instrument. Measure the absorbance at 420nm of each well.  
**Note:** The reaction time will vary depending on the level of β-galactosidase expression in the different test cultures. Absorbance values between 0.02-1.0 are within the linear range of the assay.
8. Proceed to Calculations Section.

### Microcentrifuge Tube Protocol

1. Grow cell cultures to mid-log phase (OD<sub>660</sub> of 0.5-1.0). Record the exact OD<sub>660</sub> of each culture.
2. Thaw 2X β-Galactosidase Assay Buffer on ice.
3. Add a volume of assay buffer to an equal volume of the Y-PER Reagent to make the working solution (WS). Each assay will require 350µL of the WS.
4. Prepare a blank tube containing 350µL of culture medium (no cells), 350µL of the WS and 300µL of the β-Galactosidase Assay Stop Solution.
5. Transfer 350µL of each test culture to a microcentrifuge tube.
6. Use a timer to monitor the reaction. Apply 350µL of the working reagent to each test culture and start timer.
7. Incubate the reaction tube at room temperature or 37°C (optimal) until a color change is observed.  
**Note:** The solution will become yellow within minutes depending on the amount of β-galactosidase in the sample.
8. When the yellow color appears, add 300µL of β-Galactosidase Assay Stop Solution to the reaction tube and vortex for 15 seconds. Stop the timer and record the total reaction time.
9. Remove cell debris from the reaction tube by centrifuging at 13,000 × g for 30 seconds.
10. Transfer supernatant to a cuvette and measure the absorbance at 420nm against the blank tube.  
**Note:** The reaction time will vary depending on the level of β-galactosidase expression in the test culture. Absorbance values between 0.02-1.0 are within the linear range of the assay.

### Calculations

To measure β-galactosidase activity, use the following equation:

$$\frac{1,000 \times A_{420}}{t \times V \times OD_{660}} = \beta - \text{galactosidase activity}$$

t = time (in minutes) of incubation

V = volume of cells (ml) used in the assay

### Related Thermo Scientific Products

- |              |  |
|--------------|--|
| <b>78990</b> | <b>Y-PER Yeast Protein Extraction Reagent, 500mL</b> |
| <b>75706</b> | <b>β-Galactosidase Assay Stop Solution, 25mL</b>     |
| <b>75707</b> | <b>Mammalian β-Galactosidase Assay Kit</b>           |

---

## Product Reference

1. Moriya, H. and Johnston, M. (2004). Glucose sensing and signaling in *Saccharomyces cerevisiae* through the Rgt2 glucose sensor and casein kinase I. PNAS. **101(6)**:1572-7.

This product ("Product") is warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation") and to be free from defects in material and workmanship. Unless otherwise expressly authorized in writing, Products are supplied for research use only. No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the original purchaser of the Product ("Buyer").

**No other warranties, express or implied, are granted, including without limitation, implied warranties of merchantability, fitness for any particular purpose, or non infringement. Buyer's exclusive remedy for non-conforming Products during the warranty period is limited to replacement of or refund for the non-conforming Product(s).**

There is no obligation to replace Products as the result of (i) accident, disaster or event of force majeure, (ii) misuse, fault or negligence of or by Buyer, (iii) use of the Products in a manner for which they were not designed, or (iv) improper storage and handling of the Products.

Current product instructions are available at [www.thermoscientific.com/pierce](http://www.thermoscientific.com/pierce). For a faxed copy, call 800-874-3723 or contact your local distributor.

© 2012 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.